EVOLUTION® Mechanical Dilator Sheath Set

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FM-2740-C

Cook Vascular Incorporated
Evolution® Mechanical Dilator Sheath Set

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Evolution® Mechanical Dilator Sheath Set is comprised of an inner polymer sheath connected to a handle capable of mechanically rotating the sheath and an outer telescoping polymer sheath. The inner sheath has a stainless steel tip at its distal end.

INTENDED USE

The Evolution® Mechanical Dilator Sheath Set is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects.

Use With Other Devices

The Evolution® Mechanical Dilator Sheath Set may be used in conjunction with the following catheter/lead extraction devices from COOK Vascular Incorporated:

Locking Stylet (Reference: WARNINGS);
 Byrd Stainless Steel Dilator Sheath;
 Byrd Telescoping Stainless Steel Dilator Sheath Set;
 Metal-Reinforced Flexible Dilator Sheath;
 Byrd WORK STATION™;
 NEEDLE'S EYE® Snare.

Be sure to closely follow the "Suggested Instructions for Use" for each device used.

CONTRAINDICATIONS: None known

WARNINGS

When using sheaths, do not insert sheaths over more than one lead at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

Do not use the Metal-Reinforced Flexible Dilator Sheath to apply myocardial countertraction.

When using a Locking Stylet:

Do not abandon a catheter/lead in a patient with a Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.

Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.

Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:

The item to be removed is of a dangerous shape or configuration,
 The likelihood of catheter/lead disintegration resulting in fragment embolism is high or vegetations are attached directly to the catheter/lead body.

Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.

Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.

PRECAUTIONS

Before using this device for catheter/lead removal, physicians should understand the information in:

Brodell GK, et al. Chronic transvenous pacemaker lead removal using a unique, sequential transvenous system.

Am J Cardiol 1990; 66:964-966.

Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.

Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads.

J Thorac Cardiovasc Surg 1991; 101 :989-997.

Byrd CL, et al. Intravascular lead extraction using locking wire guide stylets and sheaths.

PACE 1990; 13:1871-187 5.

Fearnot NE, et al. Intravascular lead extraction using locking wire guide stylets, sheaths, and other techniques.

PACE 1990; 13:1864-1870.

Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.

For additional copies of the preceding information, originally provided in the LEAD EXTRACTION™ System, contact COOK Vascular™ Incorporated, 800-245-4715, 724-845-8621, or FAX 724-845-2848.

Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility.

If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the non-targeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following precautions:

PREPARATIONS

Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.

Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.

Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position. Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray.

Echocardiography and cardiothoracic surgery should be rapidly available.

Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.

Establish back up pacing as necessary.

Have available an extensive collection of sheaths, Locking Stylets, stylets to unscrew active fixation leads, snares, and accessory equipment.

PROCEDURE

Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.

If using sheaths or sheath sets, including the Evolution® Mechanical Dilator Sheath Set, the following precautions should be followed:

Prior to using sheaths including the Evolution® Mechanical Dilator Sheath Set it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

The Byrd Telescoping Stainless Steel Dilator Sheath Set or Byrd Stainless Steel Dilator Sheath should only be used to minimally enter the vessel.

Be aware that insertion of a stainless steel dilator sheath over a plastic (Teflon or Polypropylene) sheath or the Evolution® Mechanical Dilator Sheath Set may damage the plastic sheath.

When advancing sheaths including the Evolution® Mechanical Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of sheaths, consider an alternate approach.

Excessive force with sheaths (including the Evolution® Mechanical Dilator Sheath Set) used intravascularly may result in damage to the vascular system requiring surgical repair.

For mechanical dissection of scar tissue using the Metal-Reinforced Flexible Dilator Sheath, turning should be limited to small, back-and-forth twisting motions (±45 degrees) to minimize mechanical stress on the device.

If the catheter/lead breaks, evaluate fragment; retrieve as indicated.

If hypotension develops, rapidly evaluate; treat as appropriate.

Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

- dislodging or damaging nontargeted catheter/lead
- chest wall hematoma
- thrombosis
- arrhythmias
- acute bacteremia
- acute hypotension
- pneumothorax
- stroke
- migrating fragment from catheter/object
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium
- cardiac tamponade
- hemothorax
- death

INSTRUCTIONS FOR USE**Suggested Instructions for Use****Evolution® Mechanical Dilator Sheath Set**

NOTE: In addition to closely following the "Suggested Instructions for Use" for this device, refer to the "Clinical Considerations for Removing Leads or Catheters: Superior Approach" shown in the next section of this insert. When using dilator sheaths or sheath sets, including the Evolution® Mechanical Dilator Sheath Set, the following precautions should be followed:

Prior to using dilator sheaths, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

The Byrd Telescoping Stainless Steel Dilator Sheath Set or Byrd Stainless Steel Dilator should only be used to minimally enter the vessel.

Be aware that insertion of a Byrd Stainless Steel Dilator Sheath over a plastic (Teflon or Polypropylene) sheath including the Evolution® Mechanical Dilator Sheath Set may damage the plastic sheath.

When advancing dilator sheaths, including the Evolution® Mechanical Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.

Excessive force with dilator sheaths, including the Evolution® Mechanical Dilator Sheath Set used intravascularly may result in damage to the vascular system requiring surgical repair.

WARNING: When using dilator sheaths or sheath sets do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

1. With the inner Evolution® Mechanical Dilator Sheath placed inside the outer sheath for telescopic action, advance the sheaths over the lead as described in Part III of "Clinical Considerations for Removing Leads or Catheters: Superior Approach" shown in the next section of this insert. To rotate the inner Evolution® Mechanical Dilator Sheath, squeeze and release the trigger activation handle. Repeat as appropriate.

2. To free the catheter/lead from tissue encapsulation refer to Part III of "Clinical Considerations for Removing Leads or Catheters: Superior Approach" shown in the next section of this insert.

Clinical Considerations for Removing Leads or Catheters:

SUPERIOR APPROACH

Clinical experience in the removal of leads from 2338 patients has identified several considerations for lead removal techniques using superior or femoral approaches. Physicians highly experienced with lead removal techniques have suggested the following considerations for removal of leads using a superior approach.

Part I. General Considerations for Accessing the Catheter/Lead

1. Surgically expose the proximal end of the indwelling catheter/lead and remove the catheter/lead from its connections (if connected).
2. Remove all suture and tie-down materials.

NOTE: Because implant procedures vary dramatically, catheters/leads may be secured in a variety of ways at the time of implant. Carefully inspect the site to ensure all suture sleeves, sutures, and tie down materials, which may be hidden or blind have been removed.

3. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.

NOTE: Do not pull on the catheter/lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to the lead may prevent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

Part II. General Considerations for Using a Locking Stylet

4. Consider passing a Locking Stylet (COOK Vascular™ Incorporated) through the inner lumen of the catheter/lead to stabilize the catheter/lead during dilation of surrounding tissues. Closely follow the "Suggested Instructions for Use" for the Locking Stylet to:

- A. Expose the inner coil of the catheter/lead
- B. Check patency of the coil lumen
- C. Determine the appropriate size of Locking Stylet based on the inner diameter of the catheter/lead
- D. Advance the Locking Stylet to the distal end of the catheter/lead
- E. Lock the Locking Stylet in place

WARNING: Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

5. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature at the proximal end of the catheter/lead, compressing the insulation against the coil and Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle or to the suture tie loop.
6. For an active fixation catheter/lead, attempt to unscrew the catheter/lead by rotating the catheter/lead (and Locking Stylet if used) counterclockwise.
7. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Locking Stylet (if used) and catheter/lead to remove it.

NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

NOTE: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

8. If the catheter/lead is not removed from the vessel with gentle pulling, dilator sheaths, including the Evolution® Mechanical Dilator Sheath Set, or other retrieval devices may help separate the catheter/lead from any tissue encapsulation as follows:

Collapse the proximal loop of the Locking Stylet (if used).

Advance dilator sheaths (or other retrieval device) over collapsed proximal loop or the proximal wire extender of the Locking Stylet (if applicable).

NOTE: If using a Wilkoff Locking Stylet, consider using a hooked stylet or small basket to aid in drawing the flattened proximal loop through the dilator sheath set (or other retrieval device).

NOTE: If using dilator sheaths or dilator sheath sets, including the Evolution® Mechanical Dilator Sheath Set (COOK Vascular™ Incorporated), consider the "General Considerations" shown below and, in addition, closely follow the "Suggested Instructions for Use" for the particular type of sheath or sheath set.

Part III. Considerations for Use of Dilator Sheaths**General Considerations**

Stainless steel, plastic (Polypropylene or Teflon), Evolution® Mechanical Dilator Sheath Sets and metal-reinforced flexible dilator sheaths are available. Stainless steel dilator sheaths or sheath sets are used to minimally enter the vessel. The stainless steel dilator sheath or sheath set is then exchanged for a plastic sheath set, Evolution® Mechanical Dilator Sheath Set (or for a metal-reinforced flexible dilator sheath, used with an outer plastic sheath). With vessel entry of the plastic Byrd Dilator Sheath Set, Evolution® Mechanical Dilator Sheath Set (or metal-reinforced flexible dilator sheath and outer plastic sheath), the sheaths are advanced over the catheter/lead in a telescoping manner, while maintaining constant tension.

Stainless steel dilator sheaths include:

Byrd Telescoping Stainless Steel Dilator Sheath Set

Byrd Stainless Steel Dilator Sheath - must be used in conjunction with the outer sheath from a plastic (Teflon or Polypropylene) sheath set

Plastic dilator sheaths include:

Byrd Dilator Sheath Set - Polypropylene

Byrd Dilator Sheath Set - Teflon

Evolution® Mechanical Dilator Sheath Set

The Metal-Reinforced Flexible Dilator Sheath must be used in conjunction with the outer sheath from a plastic (Teflon or Polypropylene) dilator sheath set.

If scar tissue makes vessel entry difficult when using a plastic Byrd Dilator Sheath Set (Polypropylene or Teflon) or Evolution® Mechanical Dilator Sheath Set, consider using the Byrd Telescoping Stainless Steel Dilator Sheath Set to obtain vessel access (refer to the "Suggested Instructions for Use" for the Byrd Telescoping Stainless Steel Dilator Sheath Set).

The size of the dilator sheaths used should be large enough that the sheaths can be advanced over the catheter/lead without causing the catheter/lead to buckle or its insulation to bunch up, but the sheaths should not be overly loose.

Always use fluoroscopic monitoring when advancing the dilator sheaths. Never use sheaths on more than one catheter or lead at a time. **Always maintain adequate tension on the catheter/lead to support the maneuvering of the dilator sheaths and properly guide them within the vessel.** With too little tension, the sheaths may rupture the vein. Too much tension may cause a myocardial avulsion.

NOTE: Follow closely the "Suggested Instructions for Use" for the particular type of sheath or sheath set being used.

- With the inner Evolution® Mechanical Dilator Sheath placed inside the appropriate outer sheath for telescopic action, insert the proximal free end of the catheter/lead into the distal end of the inner Evolution® Mechanical Dilator Sheath. Advance the catheter/lead until it completely exits the opposite (proximal) end of the sheath set.

NOTE: Consider using a hooked stylet or small basket to aid in drawing the catheter/lead through the sheath set.

- Apply adequate retracting pressure or tension on the catheter/lead and/or its wire guide or Locking Stylet (if used). **This is critical to safe passage of the sheath set over the catheter/lead. If tension is inadequate, the catheter/lead may buckle, precluding the sheath set from advancement along the appropriate path.**

NOTE: If the catheter/lead is entangled, other retrieval devices such as grasping devices, J guide wires, deflecting catheters or baskets may be helpful to grasp and straighten the catheter/lead before using dilator sheaths.

NOTE: If using stainless steel dilator sheath or sheath set, when minimal entry into the vessel has been verified, **exchange the stainless steel dilator sheath(s) for the Evolution® Mechanical Dilator Sheath Set** as described in the "Suggested Instructions for Use".

- With the catheter/lead in tension, and under fluoroscopic guidance, advance the inner sheath along the length of the catheter/lead and, if applicable, into the blood vessel. Rotate the inner Evolution® Mechanical Dilator Sheath, by squeezing and releasing the trigger activation handle. Repeat as appropriate. Advance the outer sheath over the inner sheath, keeping it within the vessel.
- In a telescopic manner, alternately advance the inner, and then outer, sheaths while maintaining adequate tension on the catheter/lead and/or its wire guide or Locking Stylet, until the catheter/lead is free of tissue restriction. In general, the inner sheath should not be advanced more than 4 cm beyond the outer sheath.

NOTE: Always monitor fluoroscopically while advancing the sheaths to avoid shearing the catheter/lead or rupturing the vessel. Continue to maintain adequate tension on the catheter/lead (or Locking Stylet). Adjust the angle of entry to keep the sheaths aligned with the catheter/lead in the vessel and to minimize sheath curvature during advancement. Rotating the inner Evolution® Mechanical Dilator Sheath using the trigger activation handle and its outer sheath during advancement may facilitate progress through exuberant scar tissue.

NOTE: If the insulation separates and exposes the inner coil, watch the tips of the outer sheath to be sure that it does not split or tear, because this could lead to vessel damage. Advance a catheter or small sheath the same size as the insulation over the bare wire, if necessary.

NOTE: If the sheaths fail to progress after initial success, or if advancing the sheaths was difficult, remove the sheaths one at a time to inspect the tips; if the tips are distorted, exchange them for a new set.

NOTE: When tracking the lead around a bend, keep the point of the tip (if angled) on the inside of the bend if applicable.

NOTE: Be careful when trying to advance the inner sheath past the ring electrode of a bipolar lead: if it gets caught, the lead may break.

- For cardiac leads, if the lead has not been freed by the time the sheaths near the myocardium, position the outer sheath so that the blunt end is directed toward the myocardium. Pull the inner sheath back several centimeters.

Apply countertraction with the outer sheath as follows:

Firmly hold the outer sheath one centimeter from the heart wall and steadily pull the Locking Stylet back; the pacing tip will be pulled into the sheath. Rotation of the sheath may help dislodge the tip.

- When the catheter/lead has become unentangled and freed from tissue attachment, remove the catheter/lead, Locking Stylet if used, and Evolution® Mechanical Dilator Sheath Set (or other retrieval device) if used.

HOW SUPPLIED

This device is supplied sterile and intended for one-time use.

REFERENCES

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. *Cleveland Clinic J Med* 1992; 59:91-92.
- Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. *Am J Cardiol* 1990;66:964-966.
- Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. *clinical cardiac Pacing*. Philadelphia: WB Saunders, 1995:491-522.
- Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.

Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. *PACE* 1990; 13:1871-1875.

Dotter CT, et al. Transluminal extraction of catheter and guide fragments from the heart and great vessels; 29 collected cases.

Am J Roentgenol 1972;111(3):467-472.

Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads. *PACE* 1993; 16:2337-2342.

Fearnot NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. *PACE* 1990;13:1864-1870.

Fisher RG, et al. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies.

Am J Roentgenol 1978;130:541-548.

Foster CJ, et al. Percutaneous removal of ventricular pacemaker electrodes using a Dormier basket. *Int J Cardiol* 1988; 21:127-134.

Smith HJ, et al. Five-years experience with intravascular lead extraction. *PACE* 1994; 17:2016-2020.

Taliercio CP, et al. Pigtail catheter for extraction of pacemaker lead (letter to the editor). *JACC* 1985;5:1020.

Trohman RG, et al. Successful percutaneous extraction of a chronic left ventricular pacing lead. *PACE* 1991; 14:1448-1451.

LEAD EXTRACTION™, and WORK STATION™ are Trademarks of COOK Vascular Incorporated. EVOLUTION®, NEEDLE'S EYE® and LIBERATOR® are Registered Trademarks of COOK Vascular Incorporated.

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Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. Cleveland Clinic J Med 1992; 59:91-92.

Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. Am J Cardiol 1990;66:964-966.

Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. clinical cardiac pacing. Philadelphia: WB Saunders, 1995:491-522.

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Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.

Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

Dotter CT, et al. Transluminal extraction of catheter and guide fragments from the heart and great vessels; 29 collected cases. Am J Roentgenol 1972;111(3):467-472.

Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads. PACE 1993; 16:2337-2342.

Fearnott NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. PACE 1990;13:1864-1870.

Fisher RG, et al. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies. Am J Roentgenol 1978;130:541-548.

Foster CJ, et al. Percutaneous removal of ventricular pacemaker electrodes using a Dormier basket. Int J Cardiol 1988; 21:127-134.

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Taliercio CP, et al. Pigtail catheter for extraction of pacemaker lead (letter to the editor). JACC 1985;5:1020.

Trohman RG, et al. Successful percutaneous extraction of a chronic left ventricular pacing lead. PACE 1991; 14:1448-1451.

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- Brodell GK, et al. Chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. Am J Cardiol 1990; 66:964-966.
 - Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
 - Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991; 101 :989-997.
 - Byrd CL, et al. Intravascular lead extraction using locking wire guide stylets and sheaths. PACE 1990; 13:1871-187 5.
 - Fearnot NE, et al. Intravascular lead extraction using locking wire guide stylets, sheaths, and other techniques. PACE 1990; 13:1864-1870.
 - Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Para obtener ejemplares adicionales de los documentos anteriores, suministrados originalmente con el sistema LEAD EXTRACTION™, póngase en contacto con COOK Vascular™ Incorporated llamando a los teléfonos 800-245-4715 o 724-845-8621, o por fax en el 724-845-2848.

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REFERENCIAS

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. Cleveland Clinic J Med 1992; 59:91-92.
- Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. Am J Cardiol 1990;66:964-966.
- Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. clinical cardiac Pacing. Philadelphia: WB Saunders, 1995:491-522.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991:101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.
- Dotter CT, et al. Transluminal extraction of catheter and guide fragments from the heart and great vessels; 29 collected cases. Am J Roentgenol 1972;111(3):467-472.
- Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads. PACE 1993; 16:2337-2342.
- Fearnott NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. PACE 1990;13:1864-1870.
- Fisher RG, et al. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies. Am J Roentgenol 1978;130:541-548.
- Foster CJ, et al. Percutaneous removal of ventricular pacemaker electrodes using a Dormier basket. Int J Cardiol 1988; 21:127-134.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Taliercio CP, et al. Pigtail catheter for extraction of pacemaker lead (letter to the editor). JACC 1985;5:1020.
- Trohman RG, et al. Successful percutaneous extraction of a chronic left ventricular pacing lead. PACE 1991; 14:1448-1451.
- LEAD EXTRACTION™ y WORK STATION™ son marcas comerciales de COOK Vascular Incorporated. EVOLUTION®, NEEDLE'S EYE® y LIBERATOR® son marcas registradas de COOK Vascular Incorporated.

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BIBLIOGRAPHIE

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. *Cleveland Clinic J Med* 1992; 59:91-92.
- Brodell GK, et al. Chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. *Am J Cardiol* 1990;66:964-966.
- Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. *Clinical Cardiac Pacing*.

Philadelphia: WB Saunders, 1995:491-522.

Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.

Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.

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Am J Roentgenol 1972;111(3):467-472.

Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads.

PACE 1993; 16:2337-2342.

Fearnout NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. *PACE* 1990;13:1864-1870.

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BIBLIOGRAFIA

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. *Cleveland Clinic J Med* 1992; 59:91-92.
- Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. *Am J Cardiol* 1990;66:964-966.
- Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. *Clinical cardiac pacing*. Philadelphia: WB Saunders, 1995:491-522.
- Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. *PACE* 1990; 13:1871-1875.
- Dotter CT, et al. Transluminal extraction of catheter and guide fragments from the heart and great vessels; 29 collected cases. *Am J Roentgenol* 1972;111(3):467-472.
- Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads. *PACE* 1993; 16:2337-2342.
- Fearnot NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. *PACE* 1990;13:1864-1870.
- Fisher RG, et al. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies. *Am J Roentgenol* 1978;130:541-548.
- Foster CJ, et al. Percutaneous removal of ventricular pacemaker electrodes using a Dormier basket. *Int J Cardiol* 1988; 21:127-134.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. *PACE* 1994; 17:2016-2020.
- Taliercio CP, et al. Pigtail catheter for extraction of pacemaker lead (letter to the editor). *JACC* 1985;5:1020.
- Trohman RG, et al. Successful percutaneous extraction of a chronic left ventricular pacing lead. *PACE* 1991; 14:1448-1451.
- LEAD EXTRACTION™ e WORK STATION™ sono marchi commerciali di COOK Vascular Incorporated. EVOLUTION®, NEEDLE'S EYE® e LIBERATOR™ sono marchi depositati di COOK Vascular Incorporated.

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- Brodell GK, et al. Chronic transvenous pacemaker lead removal using a unique, sequential transvenous system.
- Am J Cardiol 1990; 66:964-966.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads.
- J Thorac Cardiovasc Surg 1991; 101 :989-997.
- Byrd CL, et al. Intravascular lead extraction using locking wire guide stylets and sheaths.
- PACE 1990; 13:1871-187 5.
- Fearnot NE, et al. Intravascular lead extraction using locking wire guide stylets, sheaths, and other techniques. PACE 1990; 13:1864-1870.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Neem voor extra kopieën van bovenstaande informatie, oorspronkelijk meegeleverd met het LEAD EXTRACTION™ systeem, contact op met COOK Vascular™ Incorporated, 800-245-4715, 724-845-8621 of fax 724-845-2848.

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- Brodell GK, et al. Chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. Am J Cardiol 1990; 66:964-966.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991; 101 :989-997.
- Byrd CL, et al. Intravascular lead extraction using locking wire guide stylets and sheaths. PACE 1990; 13:1871-187 5.
- Fearnot NE, et al. Intravascular lead extraction using locking wire guide stylets, sheaths, and other techniques. PACE 1990; 13:1864-1870.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.

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BIBLIOGRAFIA

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. *Cleveland Clinic J Med* 1992; 59:91-92.
- Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. *Am J Cardiol* 1990;66:964-966.
- Byrd CL. Management of implant complications. In Ellenbogen KA, Kay GN, Wilkoff BL, eds. *clinical cardiac Pacing*. Philadelphia: WB Saunders, 1995:491-522.
- Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. *PACE* 1990; 13:1871-1875.
- Dotter CT, et al. Transluminal extraction of catheter and guide fragments from the heart and great vessels; 29 collected cases. *Am J Roentgenol* 1972;111(3):467-472.
- Espinosa RE, et al. The Dotter retriever and pigtail catheter: efficacy in extraction of chronic transvenous pacemaker leads. *PACE* 1993; 16:2337-2342.
- Fearnout NE, et al. Intravascular lead extraction using Locking Stylets, sheaths, and other techniques. *PACE* 1990;13:1864-1870.
- Fisher RG, et al. Evaluation of current techniques for nonsurgical removal of intravascular iatrogenic foreign bodies. *Am J Roentgenol* 1978;130:541-548.
- Foster CJ, et al. Percutaneous removal of ventricular pacemaker electrodes using a Dormier basket. *Int J Cardiol* 1988; 21:127-134.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. *PACE* 1994; 17:2016-2020.
- Taliercio CP, et al. Pigtail catheter for extraction of pacemaker lead (letter to the editor). *JACC* 1985;5:1020.
- Trohman RG, et al. Successful percutaneous extraction of a chronic left ventricular pacing lead. *PACE* 1991; 14:1448-1451.
- LEAD EXTRACTION™ e WORK STATION™ são marcas registradas da COOK Vascular Incorporated. EVOLUTION®, NEEDLE'S EYE® e LIBERATOR® são marcas registradas da COOK Vascular Incorporated.

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REFERENSER

- Brodell GK, et al. A novel approach to determining the cause of pacemaker lead failure. *Cleveland Clinic J Med* 1992; 59:91-92.
- Brodell GK, et al. chronic transvenous pacemaker lead removal using a unique, sequential transvenous system. *Am J Cardiol* 1990;66:964-966.
- Byrd CL. Management of implant complications. In: Ellenbogen KA, Kay GN, Wilkoff BL, eds. *Clinical cardiac Pacing*. Philadelphia: WB Saunders, 1995:491-522.
- Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.

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