Biliary Dilation Balloon

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Galdevejsdilatationsballon

Biliaire dilatatieballon

Ballonnet de dilatation biliaire

Gallengang-Dilatationsballon

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INTENDED USE

This device is used to dilate strictures of the biliary tree. Supplied sterile and intended for single use.

NOTES

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with balloon dilation.

Contraindications to dilation include, but are not limited to: uncooperative patient, asymptomatic strictures, inability to advance the dilation balloon through the strictured area, coagulopathy, known or suspected perforation, severe inflammation or scarring near the dilation site.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.

The duodenoscope must remain as straight as possible with the elevator in the open position when advancing or withdrawing the dilation balloon.

The entire dilation balloon must be extended beyond the tip of the duodenoscope and be fluoroscopically visualized and positioned before inflation.

The dilation balloon is used in conjunction with an inflation device or a manometer and fluid-filled syringe. **DO NOT** use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.

During withdrawal of the dilation balloon from the duodenoscope, negative pressure is mandatory to maintain balloon deflation.

A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.

Injection of contrast media during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

WARNINGS

Do not exceed the recommended balloon inflation pressure as listed on the inflation device and on the catheter tag of the dilation balloon.

Do not advance the dilation balloon if resistance is encountered. Assess the cause of resistance to determine if dilation should be re-attempted.

SYSTEM PREPARATION

1. Attach inflation device to dilation balloon.
2. Create and maintain a vacuum with inflation device. Dilation balloon is now ready for placement through accessory channel of duodenoscope.
3. Remove protective sheath from balloon. Apply a water-soluble lubricant to balloon to allow easier passage through accessory channel.

INSTRUCTIONS FOR USE (See fig.1)

I. If using Intra Ductal Exchange (IDE) port and short wire:

1. Unlock short wire from Wire Guide Locking Device and advance tip of balloon catheter over pre-positioned wire guide ensuring that wire guide exits catheter at IDE port.
2. Introduce deflated balloon into endoscope accessory channel and relock wire guide. Then advance device in short increments until endoscopically visualized exiting scope.
3. Using fluoroscopic monitoring, continue advancement into duct until desired position is achieved. **Note:** Position may be verified using fluoroscopic markers located at either end of balloon.
4. Inflate balloon to specified balloon pressure. Recommended pressures are listed on catheter tag of device.
5. Monitor balloon fluoroscopically during procedure.
6. To deflate balloon, decrease balloon pressure below 20 psi. Then create and lock in negative pressure. **Note:** Indicator on pressure gauge should be in vacuum portion of gauge during deflation.
7. Maintain negative pressure and observe balloon for complete deflation.
8. Observe balloon size. If balloon is not completely deflated, create negative pressure again as indicated in steps 6 & 7. Repeat these steps as necessary

until balloon is completely deflated. **Caution:** Balloon must be thoroughly deflated, with all fluid removed, before withdrawal. A compromised balloon may prevent removal from accessory channel, requiring removal of endoscope along with balloon.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

9. Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
10. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
11. Advance disengaged wire guide to maintain ductal access.
12. Lock wire guide into Wire Guide Locking Device.
13. To remove deflated balloon from accessory channel, straighten endoscope tip and open elevator. Maintain negative pressure while withdrawing balloon using a continuous twisting motion.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide:

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Advance tip of balloon catheter over pre-positioned long wire guide, ensuring that wire guide exits catheter at PWP.
3. Advance deflated balloon in short increments through accessory channel using standard long wire exchange technique until visualized exiting endoscope.

REFER TO STEPS 3-8 IN "SECTION I" THEN RESUME WITH STEP 4 BELOW:

4. To remove deflated balloon from accessory channel, straighten endoscope tip and open elevator. Maintain negative pressure while withdrawing balloon using a continuous twisting motion.

Upon completion of procedure, dispose of devices per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.