OASIS® One Action Stent Introduction System

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INTENDED USE

This device is intended for endoscopic biliary stent placement to drain obstructed bile ducts. This device is supplied sterile and intended for single use only.

NOTES

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that may occur with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Only Cook Cotton-Huibregtse, Cotton-Leung and ST-2 Soehendra Tannenbaum biliary stents are recommended for use with this device.

Do not use this device if the stent cannot advance through the strictured area.

SYSTEM PREPARATION

1. If applicable, pre-load desired stent and positioning sleeve onto tip of introducer.

INSTRUCTIONS FOR USE (See fig. 1)

I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire.

1. Unlock short wire from Wire Guide Locking Device and advance introducer with preloaded stent over pre-positioned wire guide ensuring that wire guide exits guiding catheter at IDE port.
2. Introduce device into endoscope accessory channel and relock wire guide.

3. With elevator open, advance device in short increments until endoscopically visualized exiting scope. **Note:** As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.
4. Advance device into appropriate duct.
5. Fluoroscopically monitor radiopaque band on tip of guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer. **Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
6. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
7. Fluoroscopically and endoscopically confirm desired stent placement.
8. Unlock wire guide from Wire Guide Locking Device.
9. Retract wire guide until it exits guiding catheter at IDE port.

Note: Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.

10. Advance disengaged wire guide to maintain ductal access.
11. Relock wire guide into Wire Guide Locking Device.
12. Retract guiding catheter into endoscope while maintaining position of pushing catheter.
13. Confirm position of stent once guiding catheter is removed.
14. Remove pushing catheter.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide.

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. **Note:** Ensure that wire guide does not exit at IDE port.

REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:

3. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
4. Confirm position of stent once guiding catheter and wire guide are removed.
5. Remove pushing catheter.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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OASIS, Cotton-Huibregtse, Cotton-Leung and Soehendra Tannenbaum are registered trademarks of Cook Endoscopy.
Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.