**COOK®**

**MEDICAL**

**Extraction Basket**

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**Panier d'extraction**

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## **INTENDED USE**

This device is used for the endoscopic removal of biliary stones and foreign bodies.

## **NOTES**

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

## **CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic procedure to be performed in conjunction with the removal of the foreign body.

If this device is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

## **POTENTIAL COMPLICATIONS**

Potential complications associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Other potential complications associated with biliary stone or foreign body removal include, but are not limited to: impaction of the object, aspiration of foreign body, localized inflammation, pressure necrosis.

## **PRECAUTIONS**

Refer to package label for the minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

If this device is to be used for the removal of biliary stones, assessment of the stone size and ampullary orifice must be made to determine the necessity of sphincterotomy. In the event sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

## **WARNING**

Device is **not** lithotripsy compatible.

When retrieving a foreign body from the upper GI tract, extreme caution must be exercised to avoid occlusion of the airway.

Surgical intervention may be required if impaction occurs.

## **INSTRUCTIONS FOR USE (See fig. 1)**

### **I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire.**

1. With basket fully retracted into sheath, unlock short wire from Wire Guide Locking Device and advance tip of basket catheter over pre-positioned wire guide ensuring that wire guide exits catheter at IDE port.
2. Insert device into endoscope accessory channel, relock wire guide into Wire Guide Locking Device, and, with elevator open, advance basket in short increments until endoscopically visualized exiting scope.
3. Endoscopically or fluoroscopically visualize extraction target and advance basket to appropriate position. Place sheath of basket beyond target before extending basket to entrap extraction target.
4. Confirm desired position of basket sheath relative to target.
5. Advance basket out of sheath by pushing forward on handle. Place basket around extraction target and slowly retract handle until entrapment is achieved. **Note:** It may be necessary to move basket back and forth along target in order to capture object.
6. Withdraw basket into sheath while maintaining entrapment.

**Caution:** If difficulty is encountered when removing basket from duct, moderate force may be applied by pulling on handle to manually fracture the stone. If passage is still restricted, surgical intervention may be necessary.

**Note:** After removing all stones, previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

7. Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
8. Retract wire guide until radiopaque distal tip of wire guide disengages from wire guide lumen.
9. Advance disengaged wire guide to maintain ductal access.
10. Lock wire guide into Wire Guide Locking Device and remove device from endoscope accessory channel.
11. Prepare or dispose of extracted specimen per institutional guidelines.

### **II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide.**

**Note:** For best results, wire guide should be kept wet.

1. With basket fully retracted into sheath, insert tip of basket catheter onto pre-positioned long wire guide, ensuring that wire guide exits catheter at PWP.
2. Advance device in short increments through accessory channel using standard long wire exchange technique until visualized exiting endoscope.

**REFER TO STEPS 3-6 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:**

3. Remove device using standard long wire exchange technique.
4. Prepare or dispose of extracted specimen per institutional guidelines.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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Soehendra is a registered trademark of Cook Medical.
Cook is a registered trademark of Cook Incorporated.
**Cook Ireland Ltd.**