**Pushing Catheter**

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**Skubbekateter**

**Pushing-katheter**

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## **INTENDED USE**

Used for endoscopic biliary stent placement to drain obstructed bile ducts. This device is supplied sterile and intended for single use only.

## **NOTES**

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

## **CONTRAINDICATIONS**

Contraindications include those specific to ERCP.

## **POTENTIAL COMPLICATIONS**

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that may occur with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

## **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Only Cook biliary and pancreatic stents are recommended for use with this device.

Do not use this device if stent cannot advance through strictured area.

## **INSTRUCTIONS FOR USE (See fig. 1)**

### **I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire:**

1. Unlock short wire from Wire Guide Locking Device and load stent and positioning sleeve over pre-positioned wire guide. **Note:** Use positioning sleeve included with this device.
2. Advance pushing catheter over pre-positioned wire guide, ensuring that wire guide exits pushing catheter at IDE port.

3. As the front flap of stent enters accessory channel, keep positioning sleeve over the back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter and wire guide, keeping it clear of accessory channel.
4. Relock wire guide.
5. With elevator open, advance device in short increments until visualized exiting endoscope.
6. Advance device into appropriate duct.
7. Fluoroscopically and endoscopically confirm desired stent placement.
8. Retract pushing catheter until IDE port is visible at Wire Guide Locking Device and resistance is felt.
9. Unlock wire guide and perform standard exchange over remaining 6 cm of pushing catheter. Once distal tip exits Wire Guide Locking Device, relock wire guide.

**II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide:**

**Note:** For best results, wire guide should be kept wet.

1. Remove stylet.
2. Load stent and positioning sleeve onto pre-positioned wire guide.
3. Advance pushing catheter over pre-positioned wire guide, ensuring that wire guide exits pushing catheter at PWP.
4. As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.

**REFER TO STEPS 5-7 IN "SECTION I", THEN RESUME WITH STEP 5 BELOW:**

5. Remove pushing catheter using standard long wire exchange technique.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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Cook is a registered trademark of Cook Incorporated.
**Cook Ireland Ltd.**