Solus™ Double Pigtail SOF-FLEX® Stent and Introduction System

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INTENDED USE

The Solus Double Pigtail SOF-FLEX Stent is used to drain obstructed biliary ducts. Each stent is recommended for use with the Solus Stent Introduction System, which includes a pushing catheter, a guiding catheter and pigtail straightener (if applicable). The Solus Stent Introduction System is included with the stent set. This device and introduction system are supplied sterile and intended for single use only.

NOTES

Do not use this device for any other purpose than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device.

If the product package is open or damaged when received, do not use the device.

Wilson-Cook devices must be stored in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass the wire guide or stent through the obstructed area.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur in conjunction with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

The coordination of the endoscope accessory channel size with compatible devices is essential in order to obtain optimal results during a procedure. Please refer to the product package label for the minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire guide device must be compatible. The Color Connection® provides instant recognition of compatible wire guides and devices through color coding on the device and package.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

CAUTIONS

A thorough understanding of the technical principles, clinical applications and risks associated with ERCP and stent placement is necessary before using this device. Solus Double Pigtail SOF-FLEX Stents and Introduction System should only be used by, or under the supervision of, physicians thoroughly trained in ERCP and stent placement.

Solus Stents must be placed under fluoroscopic monitoring.

Solus Stents must only be placed using the Solus Stent Introduction System, which is included with this stent set.

Solus Stents should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

If the wire guide or stent cannot advance through the obstructed area, do not attempt to place the stent.

The pigtail straightener is used to straighten the pigtail curls on the stent and to aid introduction of the stent into the endoscope accessory channel. The pigtail straightener is not intended for use in the accessory channel of the endoscope.

Wilson-Cook wire guides are recommended for use with this product.

Assessment must be made to determine the necessity of sphincterotomy prior to stent placement. In the event sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

WARNING

MRI compatibility of Solus Stents has not been established.

INSTRUCTIONS FOR USE

1. A complete diagnostic evaluation must be performed prior to use to determine proper stent size.
2. Upon removing the stent and its components from the package, visually inspect with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization.
3. Advance the pigtail straightener, if applicable, along the shaft of the stent to the end. (*See fig. 1*) **Note:** Care must be exercised when straightening the pigtail curls in order to avoid kinking or breaking the stent.
4. Lubricate the guiding catheter with water-soluble lubricant. Back load the stent and pigtail straightener onto the tip of the soft stent introducer. Advance the pre-loaded stent and soft stent introducer over a pre-positioned wire guide.
5. Slide the pigtail straightener, if applicable, back along the stent, straightening the second pigtail curl. Advance the stent and remove

the pigtail straightener as the stent enters the accessory channel of the endoscope. Dispose of the pigtail straightener per institutional guidelines for biohazardous medical waste. **Caution:** The pigtail straightener is not intended for use in the endoscope accessory channel.

6. Fluoroscopically monitor the radiopaque bands on the tip of the guiding catheter. The radiopaque bands on the catheter are arranged in 5 cm increments. When a sufficient length of the guiding catheter (a minimum of one radiopaque band) is above the majority of the stricture, disconnect the Luer lock fitting on the introducer. **Caution:** A sufficient portion of the guiding catheter must be above the stricture prior to stent positioning.
7. Maintain the position of the guiding catheter and advance the stent into the desired position using the pushing catheter.
8. Fluoroscopically and endoscopically confirm the desired stent placement.
9. Retract and remove both the guiding catheter and wire guide into the endoscope while maintaining the position of the pushing catheter.
10. Confirm the position of the stent once the guiding catheter and wire guide are removed.
11. Remove the pushing catheter.
12. Upon completion of the procedure, dispose of the devices per institutional guidelines for biohazardous medical waste.

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
Solus is a trademark of Wilson-Cook Medical Inc.

Color Connection is a registered trademark of Wilson-Cook Medical Inc.

Sof-Flex and Cook are registered trademarks of Cook Incorporated.

RxOnly

Wilson-Cook Medical Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

**Solus™ Double Pigtail SOF-FLEX® Stent and
Introduction System**

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18852/0105

COOK® Wilson-Cook Medical
GI Endoscopy