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INTENDED USE

The Anti-Reflux Z-Stent is a partially coated, self-expanding tubular prosthesis used to maintain patency of a malignant esophageal stricture and/or to seal a tracheoesophageal fistula across the gastroesophageal junction. The Wilson-Cook Z-Stent Introduction System is used to place the Anti-Reflux Z-Stent. This device and introduction system are supplied non-sterile and intended for single use only.

STENT DESCRIPTION

The Anti-Reflux Z-Stent is constructed of a variable number of 2 cm cages, depending on the length of the stent. When fully expanded, the stent forms a tubular shape 18 mm in diameter and is flared at both ends to allow greater anchoring within the esophageal wall. A small portion of the stent is uncoated on one end to enhance anchoring. The other end of the stent is equipped with a pressure-sensitive “windsock” that may reduce acid reflux into the esophagus and passive reflux during sleep when airways are most vulnerable and the defense mechanism against aspiration of gastric contents are not as effective.

INTRODUCTION SYSTEM DESCRIPTION

The Wilson-Cook Z-Stent Introduction System consists of a radiopaque inner catheter, an outer sheath with graduated markings, a guiding catheter, a radiopaque dilator tip, a loading funnel, a safety locking key, and a 200 cm Savary-Gilliard® wire guide.

NOTES

Do not use this device for any purpose other than the stated intended use.

If the product package is open or damaged when received, do not use this device.

Wilson-Cook devices must be stored in a dry location, away from temperature extremes.

The Anti-Reflux Z-Stent may not be sterilized by any method or soaked in any solution.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to the placement of the Anti-Reflux Z-Stent include, but are not limited to: total esophageal obstruction, strictures that cannot be dilated to 10 mm, placement requiring positioning of the Z-Stent within 2 cm of the cricopharynx, tumors located proximal to mid-esophagus, surgical resection candidates, patients with a perforated esophagus, patients with a hiatal hernia, placement in actively bleeding tumors, benign disease, uncooperative patient, coagulopathy, tracheal compression, recent myocardial infarction, cervical arthritis with fixed cervical spine, large tumor

mass occupying the mediastinum, non-obstructing tumor, gastric outlet obstruction, necrotic esophageal mucosa and acutely angled stenosis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper GI endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur with the placement of an esophageal prosthesis include, but are not limited to: stent misplacement; stent migration; tumor ingrowth/overgrowth; esophageal ulceration; esophageal erosion; nausea; chest or retrosternal pain; foreign body sensation; food bolus impaction; gas-bloat; sensitivity to metal component(s); fistula involving trachea, bronchi, or pleural space; intestinal obstruction secondary to migration; mediastinitis/peritonitis; airway compression/tracheal obstruction.

PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine proper Z-Stent size.

Prior to loading the Z-Stent, ensure that the loading string is not twisted.

CAUTIONS

A thorough understanding of the principles, clinical applications and risks associated with upper GI endoscopy and esophageal stent placement is necessary before using this device. The Anti-Reflux Z-Stent should only be placed by, or under the supervision of, physicians thoroughly trained in the procedure.

The Anti-Reflux Z-Stent should be placed using fluoroscopic monitoring.

The Anti-Reflux Z-Stent should only be placed with the Wilson-Cook Z-Stent Introduction System, which is provided with each stent.

Prior to advancing the Z-Stent Introduction System into the esophagus, the area to be stented should be dilated to a minimum of 10 mm and a maximum of 14 mm. If the area is dilated greater than 14 mm, the stent may migrate.

The Anti-Reflux Z-Stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After Anti-Reflux Z-Stent placement, patients should be instructed to chew their food well, or to eat soft or pureed food.

The use of rigid type dilator for post-placement dilation is not recommended, as this may dislodge the stent.

After Anti-Reflux Z-Stent placement, alternate methods of treatment, such as chemotherapy and irradiation, **should not** be administered as this may

increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

A minimum stent length of 4 cm longer than the lesion is recommended to bridge the stricture with 2-3 cm (1-2 cages) above and 2-3 cm (1-2 cages) below the affected area. When placing the Anti-Reflux Z-Stent across the gastroesophageal junction, the minimum stent length should be 6 cm greater than the length of the lesion, allowing placement of 4 cm above and 2 cm below the stricture.

Long-term patency of the Anti-Reflux Z-Stent has not been established. Periodic evaluation of the device is advised.

WARNINGS

The Anti-Reflux Z-Stent is **not intended to be removed**. Attempts to remove the stent after placement may cause damage to the esophageal mucosa. The Anti-Reflux Z-Stent is considered to be a permanent implant.

MRI compatibility of the Anti-Reflux Z-Stent has not been established.

EQUIPMENT REQUIRED

- Anti-Reflux Z-Stent with Introduction System
- Quantum TTC® Esophageal Balloon Dilators
- or
- Savary-Gilliard® Dilators

PATIENT PREPARATION

1. A complete diagnostic evaluation should be performed prior to placement to determine the proper Z-Stent length.
2. Upon removing the Z-Stent and introduction components from the packages, visually inspect the devices with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization.
3. Dilate the stricture to a minimum of 10 mm and a maximum of 14 mm using Quantum TTC Balloon Dilators or Savary-Gilliard Dilators.
4. Endoscopically locate and mark the upper and lower margins of the lesion with a radiopaque marker. **Caution:** Accurate marking of the lesion borders is essential for proper Z-Stent placement. **Note:** Water-soluble markers will dissipate quickly.
5. Place the Savary-Gilliard wire guide, floppy tip first, through the endoscope accessory channel until the tip is fluoroscopically visualized in the stomach. Leave the wire guide in place and remove the endoscope.

SYSTEM PREPARATION

1. Ensure the guiding catheter fitting is completely snapped into the fitting at the end of the inner pusher.

2. Apply the enclosed tube of water-soluble lubricant to the inside of the loading funnel flare and the inner rim of the outer sheath. **Note:** Prior to loading the Z-Stent, ensure that the loading string is not twisted.
3. Load the stent into the introducer by grasping the funnel and stent in one hand and pulling the thumb-ring attached to the suture with the other hand. Continue pulling the thumb-ring until the stent is fully loaded into the introducer. **Note:** The stent should be guided into the funnel, if necessary, to facilitate the loading into the introducer.
4. After the stent is fully loaded into the outer sheath, disengage the loading funnel from the outer sheath and remove the loading funnel from the guiding catheter.
5. Gently push into the outer sheath any part of the stent that is not yet fully loaded. Pull the guiding catheter fitting to disengage it from the inner catheter. Advance the dilator tip into the outer sheath. The dilator tip should fit securely within the outer sheath. The tip may advance the stent approximately 2 cm into the outer sheath.
6. Place the proximal end of the previously positioned Savary-Gilliard wire guide through the guiding catheter. Gently advance the assembled system over the wire guide until the wire guide exits the fitting on the guiding catheter.
7. Remove the locking key and advance the inner pusher to advance the stent through the overtube until the anti-reflux sleeve is fully compressed against the dilator tip. **Note:** It may be necessary to hold the dilator tip to prevent it from being deployed while advancing the stent.

INSTRUCTIONS FOR PLACEMENT

1. Using fluoroscopic guidance, advance the assembled introduction system over the pre-positioned Savary-Gilliard wire guide. Using the collapsed cages as reference points, fluoroscopically align the stent within the marked lesion. If placed in the body of the esophagus, a stent length allowing 2-3 cm (1-2 cages) above and 2-3 cm (1-2 cages) below the stricture is appropriate. If placed across the GE junction, the additional length allows 4 cm (2 cages) above and 2 cm (1 cage) below the stricture.
2. Fluoroscopically verify that the stent position within the lesion has been maintained. Reposition the introducer system as necessary.
3. Using continuous fluoroscopic monitoring, **SLOWLY** deploy the **first** cage of the stent by holding the inner catheter in place while **slowly** withdrawing the outer sheath. **Caution:** Failure to deploy the stent **slowly** may cause improper positioning. **Note:** During deployment of the first cage, the stent position may be adjusted by **slowly** advancing or withdrawing the entire introduction system. Continue deployment,

one cage at a time, by **slowly** withdrawing the outer sheath, until stent is completely deployed.

4. Cut the suture attached to the thumb-ring and remove the thumb-ring, leaving the suture in place through the stent.
5. Withdraw the introduction system, leaving the loading suture in place through the stent. **Caution:** If the loading suture becomes compressed against the guiding catheter and outer sheath during removal of the introduction system, reengage the guiding catheter into the inner catheter of the sheath to free the suture. Do **not** remove the introduction system while the suture is caught, as this may cause inadvertent displacement of the stent.
6. Introduce the endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent. **Caution:** Do not introduce the endoscope into the stent as displacement may occur. **Caution:** Tumor impingement may prevent the stent from immediately achieving the maximum 18 mm diameter. **Note:** The loading suture may be used as an aid in stent repositioning, if necessary.
7. After proper stent position is verified endoscopically, continue to endoscopically monitor the stent while gently pulling on one end of the loading suture until the suture is completely removed.
8. Upon completion of the procedure, dispose of the introduction system components per institutional guidelines for biohazardous medical waste.

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Z-Stent, Savary-Gilliard and Quantum TTC are registered trademarks of Wilson-Cook Medical Inc.

Cook is a registered trademark of Cook Incorporated.

Wilson-Cook Medical Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

**Esophageal Z-Stent[®]
With Dua Anti-Reflux Valve**

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18869/0105

COOK[®] Wilson-Cook Medical
GI Endoscopy