COOK®

MEDICAL

Memory™ Polyp Retrieval Baskets

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INTENDED USE

This device is used for the endoscopic removal of polyps after polypectomy.

NOTES

This device is supplied non-sterile and reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in conjunction with removal of polyps.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket removal of a polyp include, but are not limited to: localized inflammation.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

INSTRUCTIONS FOR USE

1. With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope.
2. Endoscopically visualize polyp to be removed and advance basket to appropriate position. Place sheath of basket beyond target, then extend basket to entrap polyp. **Note:** It may be necessary to move basket back and forth along target in order to capture polyp.
3. Withdraw basket into sheath while maintaining entrapment. Support polyp against endoscope tip, and withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of polyp. While withdrawing device from endoscope, wipe secretions from outer sheath.
4. Prepare extracted specimen per institutional guidelines.
5. Upon completion of procedure, clean and prepare device for future use following cleaning instructions. During cleaning, inspect integrity and function

of device to determine advisability of reuse. **If kinks, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.**

CLEANING INSTRUCTIONS

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaning solution to Luer lock fitting on device, then irrigate.
2. Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.
4. Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.
5. Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
6. Remove device from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended basket thoroughly with clean, running water.
7. Ultrasonically clean device, with basket extended, for minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove device from ultrasonic cleaner, then thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.
9. Gently wipe external portions of device dry with a lint-free cloth. Air dry lumen and basket by hanging device, with basket extended, or use medical grade (oil-free) forced air through Luer lock fitting on handle.
10. Package device for sterilization following AAMI recommended practices.
11. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER

Prevacuum Pressure
Relative Humidity
Humidity Dwell
Chamber Temperature
Exposure Pressure
Exposure Dwell
ETO Concentration
Post Vacuums
Post Vacuum Pressure
Air Wash²
Aeration Temperature
Aeration Dwell
Air Changes During Aeration

SET POINT

2.0 psia (25.9" Hg Vac)
60% RH
20 minutes
54.4°C (130°F)
To Be Determined¹
2 hours
600 mg/L
3
2.0 psia (25.9" Hg Vac)
5 minutes
54.4°C (130°F)
12 hours
90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space.

AUTOCLAVE CYCLE PARAMETERS

PARAMETER

Cycle Type
Prevacuum Pressure
Chamber Exposure Temperature
Chamber Exposure Pressure
Exposure Dwell

SET POINT

Prevacuum
2.5 psia (24.8" Hg Vac)
132.2°C (270°F)
41.8 psia (27.1 psig)
5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No. 12 - 1994.

"Sterilization, Part 1 - Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.


Memory is a trademark of Cook Medical.

Cook is a registered trademark of Cook Incorporated.

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

RxOnly

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

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