Howell Biliary Aspiration Needle
Agulha de Aspiração Biliar Howell
Aguja para Aspiración Biliar Howell

18905/0105

COOK[®] Wilson-Cook Medical
GI Endoscopy

English	1
Português/Portuguese	3
Español/Spanish	6

INTENDED USE

The Howell Biliary Aspiration Needle is used for aspiration biopsy in the biliary system. This device is supplied sterile and intended for single use only.

NOTES

Do not use this device for any purpose other than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device.

If the product package is opened or damaged when received, do not use the device.

Wilson-Cook devices must be stored in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTION

The coordination of the endoscope accessory channel size with compatible devices is essential in order to obtain optimal results during a procedure. The Howell Biliary Aspiration Needle catheter is 7.0 French and requires a minimum channel size of 2.8 mm.

CAUTIONS

A thorough understanding of the technical principles, clinical applications and risks associated with ERCP and aspiration biopsy is necessary before using this device. The Howell Biliary Aspiration Needle should only be used by, or under the supervision of, physicians thoroughly trained in therapeutic ERCP.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

The needle must be fully retracted into the sheath prior to introduction into, advancement through or withdrawal from the endoscope. Failure to retract the needle may result in damage to the endoscope.

INSTRUCTIONS FOR USE

1. Upon removing the device from the package, visually inspect with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization.
 2. Uncoil and straighten the catheter completely.
 3. Prior to inserting the device into the endoscope, retract the needle completely into the sheath. Ensure the needle fully retracts into the sheath and that no sharp tip extends beyond the ball tip.
Caution: Failure to retract the needle prior to introduction into, advancement through, and withdrawal from the endoscope may result in damage to the endoscope.
 4. Introduce the device into the accessory channel of the endoscope and advance in short increments to the desired biopsy site. Do not straighten the pre-bend more than is necessary for insertion down the channel. **Note:** A large bore endoscope is preferred to minimize straightening.
 5. Position the endoscope tip away from and slightly above the ampulla.
 6. Place the ball tip, with the needle retracted, on the lower lip of the biliary orifice. Advance the catheter into the cannulating position. Proceed with cannulation of the biliary system.
 7. Once cannulation is achieved, advance the ball tip to the appropriate position for aspiration biopsy.
 8. Carefully advance the needle out of the sheath. With the needle in the desired location, remove the stylet from the Luer lock fitting on the handle.
 9. Attach a 20 cc syringe to the Luer lock fitting, then aspirate. Continue aspiration and gently move the needle in small increments back and forth within the biopsy site. **Note:** Several biopsy locations may be chosen and aspiration repeated without withdrawing the needle from the bile duct. Between each placement of the needle for biopsy, suction should be released.
 10. When aspiration is complete, withdraw the needle from the biopsy site and into the ball tip by pulling back on the Luer slip fitting, then withdraw the device from the endoscope.
 11. Prepare the specimen for examination per institutional guidelines.
 12. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.
- 

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© 2005 Wilson-Cook Medical Inc.
Cook is a registered trademark of Cook Incorporated.
Wilson-Cook Medical Inc.
Cook Ireland Ltd.