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18915/0205

COOK[®] Wilson-Cook Medical
GI Endoscopy

English	1
Português/Portuguese	6
Español/Spanish	11

INTENDED USE

The Wilson-Cook Esophageal Z-Stent, Fully Coated or with Uncoated Flanges, is a self-expanding tubular prosthesis used to maintain patency of a malignant esophageal stricture and/or to seal a tracheoesophageal fistula. The Wilson-Cook Z Speed Introduction System is used to place the Esophageal Z-Stent, Fully Coated or with Uncoated Flanges. This device and introduction system are supplied non-sterile and intended for single use only.

STENT DESCRIPTION

The Wilson-Cook Esophageal Z-Stent, Fully Coated or with Uncoated Flanges, is constructed of a series of cages, depending on the length of the stent. When fully expanded, the stent forms a tubular shape that is flared at both ends to allow greater anchoring within the esophageal wall. Refer to product label for specific stent lengths and diameters.

INTRODUCTION SYSTEM DESCRIPTION

The Wilson-Cook Z Speed Introduction System consists of a radiopaque dilator, dilator locking adapter, internal sheath, and handle with flush port. The stent loading cartridge is preloaded with an Esophageal Z-Stent, Fully Coated or with Uncoated Flanges, and attached to a clear external sheath with positioning reference markers. (*See fig. 1*)

- | | |
|--------------------------------------|----------------------------|
| a) Z-Stent | g) Dilator locking adapter |
| b) Stent loading cartridge | h) Release key |
| c) Proximal cartridge fitting | i) Dilator tip |
| d) Steering cap with release sutures | j) Handle |
| e) Clear external sheath | k) Flush port |
| f) Internal sheath | |

NOTES

Do not use this device for any purpose other than the stated intended use.

If the product package is open or damaged when received, do not use this device.

Wilson-Cook devices must be stored in a dry location, away from temperature extremes.

The Esophageal Z-Stent may **not** be sterilized by any method or soaked in any solution.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to the placement of the Esophageal Z-Stent include,

but are not limited to: total esophageal obstruction, strictures that cannot be dilated to a minimum of 10 mm, placement requiring positioning of the Esophageal Z-Stent within 2 cm of the cricopharynx, surgical resection candidates, patients with a perforated esophagus, placement in actively bleeding tumors, benign disease.

Relative contraindications include, but are not limited to: uncooperative patient, coagulopathy, tracheal compression, recent myocardial infarction, cervical arthritis with fixed cervical spine, large tumor mass occupying the mediastinum, non-obstructing tumor, gastric outlet obstruction, necrotic esophageal mucosa, and acutely angled stenosis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper GI endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, reflux, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur with the placement of an esophageal prosthesis include, but are not limited to: stent misplacement; stent migration; tumor ingrowth/overgrowth; esophageal ulceration; esophageal erosion; nausea; chest or retrosternal pain; foreign body sensation; food bolus impaction; gas-bloat; sensitivity to metal component(s); fistula involving trachea, bronchi, or pleural space; intestinal obstruction secondary to migration; mediastinitis/peritonitis; airway compression/tracheal obstruction.

PRECAUTIONS

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

Prior to loading the Esophageal Z-Stent, ensure that the loading string is not twisted.

CAUTIONS

A thorough understanding of the technical principles, clinical applications and risks associated with upper GI endoscopy and esophageal stent placement is necessary before using this device. The Esophageal Z-Stent should only be placed by, or under the supervision of, physicians thoroughly trained in the procedure.

The Esophageal Z-Stent should be placed using fluoroscopic monitoring.

The Esophageal Z-Stent must only be placed with the Wilson-Cook Z Speed Introduction System, which is provided with each stent.

Prior to advancing the Esophageal Z-Stent Introduction System into the esophagus, the area to be stented should be dilated to a minimum of 10 mm and a maximum of 14 mm. If the area is dilated greater than 14 mm, the stent may migrate.

The Esophageal Z-Stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After Esophageal Z-Stent placement, alternative methods of treatment such as chemotherapy and irradiation **should not** be administered as this may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

A minimum stent length of 4 cm longer than the lesion is recommended to bridge the stricture with 2-3 cm (1-2 cages) above and 2-3 cm (1-2 cages) below the affected area. When placing an Esophageal Z-Stent across the gastroesophageal junction, the minimum stent length should be 6 cm greater than the length of the lesion, allowing placement of 4 cm above and 2 cm below the stricture.

Long-term patency of the Esophageal Z-Stent has not been established. Periodic evaluation of the device is advised.

After Esophageal Z-Stent placement, patients should be instructed to chew their food well, or to eat soft or pureed food.

The use of rigid type dilator for post-placement dilation is not recommended, as this may dislodge the stent.

WARNINGS

The Esophageal Z-Stent is **not intended to be removed**. Attempts to remove the stent after placement may cause damage to the esophageal mucosa. The Esophageal Z-Stent is considered to be a permanent implant.

MRI compatibility of the Esophageal Z-Stent has not been established.

EQUIPMENT REQUIRED

- Esophageal Z-Stent, Fully Coated or with Uncoated Flanges, with Z Speed Introduction System
- Wilson-Cook Esophageal Balloon Dilators
or
- Savary-Gilliard® Dilators
- 5.5 French catheter (ERCP-1 recommended)
- .035 inch wire guide (THSF recommended)

PATIENT PREPARATION

1. A complete diagnostic evaluation should be performed prior to placement to determine the proper stent size.
2. Upon removing the Esophageal Z-Stent and introduction components from the packages, visually inspect the devices with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization.

3. Dilate the stricture to a minimum of 10 mm and a maximum of 14 mm using a Wilson-Cook Balloon Dilator or a Savary-Gilliard Dilator.
4. Endoscopically locate and mark the upper and lower margins of the lesion with a radiopaque marker. **Caution:** Accurate marking of the lesion borders is essential for proper Z-Stent placement.
Note: Water-soluble markers will dissipate quickly.
5. Introduce a 5.5 French catheter through the endoscope accessory channel until it is fluoroscopically visualized in the stomach.
6. Place .035 inch wire guide, floppy tip first, through the catheter until it is fluoroscopically visualized in the stomach.
7. Leave the wire guide in place and remove the catheter and the endoscope.

SYSTEM PREPARATION

1. Remove protective caps from proximal cartridge fitting and clear external sheath. Remove shaping wire from dilator tip.
2. Introduce the dilator tip (C) into the stent loading cartridge through the stent until the dilator tip enters the steering cap and the proximal end of the external sheath.
3. Push the dilator, steering cap and stent through the external sheath until the dilator locking adapter reaches the proximal cartridge fitting and lock them together.
4. Cut one of the release sutures and remove the steering cap and the suture. (*See fig. 2*)
5. Lubricate distal dilator tip and clear external sheath.
6. Loading has now been completed and the stent is ready for placement.
Note: Do not use the product if there is any doubt as to whether the stent has been loaded properly.

INSTRUCTIONS FOR PLACEMENT

1. Using fluoroscopic guidance, introduce and advance the introduction system over the pre-positioned wire guide and into the esophagus.
Note: The stent is visible under fluoroscopy. (*See fig. 3*)
2. Using the position markers on the exterior sheath as a reference, advance the introduction system through the stricture. Position the distal portion of the stent 2 to 3 cm or 1-2 cages past the stricture. The stent must cover the tumor completely. If this is not possible, add a proximal or distal stent extension. **Note:** If the tumor area is not located satisfactorily or there is any doubt where to deploy the stent, place catheter along the introducer system. Inject contrast medium to visualize the proximal end of the obstruction.

3. With the distal end of the stent at desired position, remove the release key from the dilator locking adapter.
4. To deploy the stent, hold the handle and internal sheath stationary while pulling the loading cartridge and external sheath towards the handle and internal sheath. (*See fig. 4*)
5. The stent has been fully deployed when the loading cartridge reaches the internal handle. **Note:** After full deployment of the stent, further stent maneuvering is not possible. (*See fig. 5*)
6. Introduce the endoscope and advance to the upper margin of the stent to endoscopically confirm position and patency of the stent. **Caution:** Do not introduce the endoscope into the stent as displacement may occur. **Caution:** Tumor impingement may prevent the stent from immediately achieving its maximum diameter. A dilating balloon may be used to expand the stent to its maximum diameter.
7. Upon completion of the procedure, remove the stent introducer sheath, dilator assembly and the wire guide and dispose of per institutional guidelines for biohazardous medical waste.

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



Z Speed is a trademark of Wilson-Cook Medical Inc.

Z-Stent and Savary-Gilliard are registered trademarks of Wilson-Cook Medical Inc.

Cook is a registered trademark of Cook Incorporated.

Wilson-Cook Medical Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

**Esophageal Z-Stent[®] Fully Coated or
Uncoated Flanges Endoprosthesis System
with Z Speed[™] Introducer**

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