**Soehendra[®] Stent Retriever /
Soehendra Universal Catheter**

**Soehendra[®] Extractor de Stents /
Soehendra Cateter Universal**

**Extractor de Stents Soehendra[®] /
Catéter Universal Soehendra**

18916/0205

COOK[®] Wilson-Cook Medical
GI Endoscopy

English	1
Português/Portuguese	5
Español/Spanish	10

INTENDED USE

The Soehendra Universal Catheter and Soehendra Stent Retriever are used in conjunction to cannulate and remove stents from the biliary and pancreatic ducts while maintaining wire guide placement. The Soehendra Universal Catheter is supplied sterile and is reusable if the device integrity is intact. The Soehendra Stent Retriever is supplied sterile and intended for single use only.

NOTES

Do not use this device for any purpose other than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device.

If the product package is open or damaged when received, do not use this device.

Wilson-Cook devices must be stored in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

The coordination of endoscope accessory channel size with compatible devices is essential in obtaining optimal results during a procedure. The Soehendra Universal Catheter is 6.0 FR and requires a minimum channel size of 2.8 mm. The Soehendra Stent Retriever is available in various French sizes depending upon the size stent to be retrieved. The following chart provides rapid reference for assuring compatibility of the channel and stent with the desired retriever.

STENT TO BE RETRIEVED	ENDOSCOPE CHANNEL
French Size	Minimum Size
5.0	2.8 mm
7.0	2.8 mm
8.5	3.2 mm
10.0	3.7 mm
11.5	4.2 mm

Wilson-Cook wire guides are recommended for use with this device.

Injection of contrast during ERCP must be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

The reusability of a device depends in large part on the care of the device by the user. Factors involved in prolonging the life of the Soehendra Universal Catheter include, but are not limited to: advancing through the endoscope accessory channel in short increments, gently withdrawing the device from the channel, avoiding loops and/or kinks in the catheter, thorough cleaning following the instructions included in this booklet.

INSTRUCTIONS FOR USE

1. Upon removing the Soehendra Universal Catheter from the package and prior to each use, visually inspect with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization if an abnormality is noted upon initial inspection of the device.
2. Upon removing the Soehendra Stent Retriever from the package, visually inspect with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization.
Note: If the product package has been opened inadvertently and the device not used, the device may be resterilized **one time only**. Refer to the table at the end of these instructions for the recommended ETO sterilization cycle parameters. Package the device for sterilization following AAMI recommended practices.
3. Introduce and advance the Soehendra Universal Catheter in small increments through the endoscope accessory channel. **Note:** Optional when using the Soehendra Stent Retriever with extended tip.
4. Advance the appropriate wire guide through the Soehendra Universal Catheter in small increments until the wire guide is completely through the stent and well into the duct.
5. Remove the catheter, then clean and prepare it for future use following the cleaning instructions that follow. During cleaning, inspect the integrity and function of the device to determine advisability of reuse. If kinks, bends or breaks exist, do not use the device. Dispose of the device per institutional guidelines for biohazardous medical waste.
6. Introduce and advance the threaded end of the Soehendra Stent Retriever over the pre-positioned wire guide until it reaches the stent. (See fig. 1) **Note:** Soehendra Stent Retriever with extended tip may be used prior to wire guide placement.
7. Keeping the end of the Soehendra Stent Retriever that is outside of the scope straight, turn the stent retriever handle clockwise until the threaded end completely connects to the stent. (See fig. 2) **Note:** Maximum thread engagement is recommended.

8. Once the stent is securely connected to the retriever, open the elevator of the endoscope and slowly pull the stent and Soehendra Stent Retriever out of the channel making sure the wire guide is left in place. **Note:** Dislodgement of the wire guide from the duct can be prevented by simultaneously pushing the wire guide and pulling the stent and stent retriever out of the endoscope. **Caution:** If the stent detaches from the Soehendra Stent Retriever while in the channel, push the stent down the scope and into the elevator with the stent retriever. Re-thread the stent onto the retriever, then slowly pull both out of the endoscope.
9. Upon completion of the procedure, dispose of the Soehendra Stent Retriever per institutional guidelines for biohazardous medical waste.

UNIVERSAL CATHETER CLEANING INSTRUCTIONS

Caution: When handling, coil the Soehendra Universal Catheter in a minimum of 8-inch diameter to prevent product damage.

1. While withdrawing the Soehendra Universal Catheter from the endoscope, wipe the excess secretions from the outer sheath.
2. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to the Luer lock fitting on the catheter, then irrigate the inner lumen of the catheter.
3. Submerge the solution-filled catheter in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
4. Beginning at the Luer lock fitting and working toward the tip, wipe the catheter with a textured cloth soaked in enzymatic solution.
5. Irrigate the inner lumen of the catheter by attaching a syringe containing a minimum of 20 cc of enzymatic solution to the Luer lock fitting on the catheter. Repeat this step until all visible secretions are removed.
6. Remove the catheter from the enzymatic solution, then thoroughly rinse the inner lumen of the catheter using a syringe filled with clean water. Rinse the outer portions of the catheter thoroughly with clean, running water.
7. Ultrasonically clean the catheter for a minimum of 15 minutes at 25°C (77°F) or per the ultrasonic cleaner manufacturer's instructions.
8. Remove the catheter from the ultrasonic cleaner, then thoroughly rinse the inner lumen with a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse the outer portions of the catheter with clean, running water.
9. Gently wipe the external portions of the catheter dry with a lint-free cloth. Hang the device to air dry the inner lumen, or use medical grade (oil free) forced air through the Luer lock fitting.
10. Package the device for sterilization following AAMI recommended practices. **Note:** AAMI stands for Association for the Advancement of Medical Instrumentation.

11. Sterilize the device using the set points for the desired sterilization method shown in the following tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (-25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (-25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

AUTOCLAVE CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (-24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No. 12 - 1994.

"Sterilization, Part 1 - Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

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

“Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers”, *AAMI Technical Information Report*, No. 12 - 1994.

“Sterilization, Part 1 - Good Hospital Practices”, *AAMI Standards and Recommended Practices*, Volumen 1.1 - 1995.

Soehendra is a registered trademark of Wilson-Cook Medical Inc.

Cook is a registered trademark of Cook Incorporated.

**Soehendra Stent
Retriever Only**

Wilson-Cook Medical Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

RxOnly

STERILE	EO
----------------	-----------

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

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