**COOK®**

**MEDICAL**

**Classic ERCP Catheters**

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## INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

## NOTES

This device is reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

## CONTRAINDICATIONS

Those specific to ERCP.

## POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Reusability of device depends in large part on care of device by user. Factors involved in prolonging life of this device include, but are not limited to: advancing through endoscope accessory channel in short increments, gently withdrawing device from channel, avoiding loops and/or kinks in catheter, coiling device in a minimum 8-inch (20 cm) diameter, thorough cleaning following instructions included in this booklet.

## INSTRUCTIONS FOR USE

1. Gently advance stylet completely into catheter. Lock stylet cap into catheter fitting.
2. Introduce catheter into endoscope and advance in short increments until device is endoscopically visualized exiting scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm position of device by attaching a syringe to side port of catheter fitting.
4. To insert a compatible wire guide, remove stylet from catheter and advance wire guide prior to removing catheter from endoscope. **Note:** Before each introduction of a wire guide into a device, wet wire guide in a sterile water or saline bath. A contrast-filled catheter may make wire guide difficult to advance. Flush contrast from catheter with sterile water or saline before each introduction of wire guide, observing necessary cautions if device is positioned in pancreatic duct.

5. Upon completion of procedure, clean and prepare device for future use following cleaning instructions that follow. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use and dispose of device per institutional guidelines for biohazardous medical waste.

## **CLEANING INSTRUCTIONS**

1. While withdrawing catheter from endoscope, wipe excess secretions from outer sheath.
2. Immediately after use, remove stylet from catheter and attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on catheter. Block open port on catheter fitting, then irrigate inner lumen of catheter.
3. Submerge solution-filled catheter and stylet in solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
4. Beginning at Luer lock fitting and working toward tip, wipe device with a textured cloth soaked in enzymatic solution. Beginning at hub and working downward, wipe stylet.
5. Irrigate inner lumen of catheter by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on catheter. Repeat until all visible secretions are removed.
6. Remove catheter and stylet from enzymatic solution, then thoroughly rinse inner lumen of catheter using a syringe filled with clean water. Rinse stylet and outer portions of catheter thoroughly with clean, running water.
7. Ultrasonically clean catheter and stylet for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove catheter and stylet from ultrasonic cleaner, then thoroughly rinse inner lumen of catheter with a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse stylet and outer portions of catheter with clean, running water.
9. Gently wipe external portions of stylet and catheter dry with a lint-free cloth. Air dry inner lumen of catheter by hanging device for three hours, or by using medical grade (oil free) forced air through Luer lock fitting.
10. Insert stylet into catheter, leaving approximately 1 cm extending from Luer lock fitting.
11. Package device for sterilization following AAMI recommended practices.
12. Sterilize device using set points for desired sterilization method shown in tables.

## ETO CYCLE PARAMETERS

### PARAMETER

Prevacuum Pressure  
Relative Humidity  
Humidity Dwell  
Chamber Temperature  
Exposure Pressure  
Exposure Dwell  
ETO Concentration  
Post Vacuums  
Post Vacuum Pressure  
Air Wash<sup>2</sup>  
Aeration Temperature  
Aeration Dwell  
Air Changes During Aeration

### SET POINT

2.0 psia (25.9" Hg Vac)  
60% RH  
20 minutes  
54.4°C (130°F)  
To Be Determined<sup>1</sup>  
2 hours  
600 mg/L  
3  
2.0 psia (25.9" Hg Vac)  
5 minutes  
54.4°C (130°F)  
12 hours  
90 per hour

<sup>1</sup> To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

<sup>2</sup> Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

**Note:** ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

## AUTOCLAVE CYCLE PARAMETERS

### PARAMETER

Cycle Type  
Prevacuum Pressure  
Chamber Exposure Temperature  
Chamber Exposure Pressure  
Exposure Dwell

### SET POINT

Prevacuum  
2.5 psia (24.8" Hg Vac)  
132.2°C (270°F)  
41.8 psia (27.1 psig)  
5 minutes

**Note:** Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No. 12 - 1994.


"Sterilization, Part 1 - Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

Cook is a registered trademark of Cook Incorporated.

**Wilson-Cook Medical, Inc.**  
4900 Bethania Station Road  
Winston-Salem, North Carolina 27105  
USA

**RxOnly**

**Cook Ireland Ltd.**  
O'Halloran Road  
National Technology Park  
Limerick  
Ireland

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