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## **INTENDED USE**

This device is used in conjunction with a single or dual plug bipolar electro-surgical generator to endoscopically provide hemostasis throughout the gastrointestinal tract.

## **NOTES**

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

## **CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired coagulation site.

## **POTENTIAL COMPLICATIONS**

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur with the use of the bipolar hemostatic probe include, but are not limited to: transmural burns, thermal injury to the patient, explosion.

## **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

This device meets the applicable portions of the standard governing electro-surgical devices (ANSI/AAMI HF18).

Before using this device, follow recommendations provided by electro-surgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode. Ensure a proper path from patient return electrode to electro-surgical unit is maintained throughout procedure.

Bipolar coagulation should be performed under direct endoscopic visualization.

For use with bipolar generators not to exceed 1.4 kVp-p.

## **WARNINGS**

This device should never be used in a monopolar mode of operation. Do not plug into the patient grounding pad receptacle.

Do not pass probe through an endoscope whose distal end is deflected at an angle more than approximately 15 degree.

There is the possibility of explosion if probe is used in the presence of flammable anesthetics or naturally produced methane gas.

### GENERATOR COMPATIBILITY

Order Number	Generator
BCP-7A	Valley Lab / Everest Medical BiCOAG
BCP-10A	Valley Lab / Everest Medical BiCOAG
BCP-7B	BiCap and Bicap II / Bard / Microvasive
BCP-10B	BiCap and Bicap II / Bard / Microvasive

### INSTRUCTIONS FOR USE

1. With electrosurgical unit off, prepare equipment.
2. Connect flushing port Luer connector to a distilled water source (syringe or pump) to facilitate irrigation of site that is to be coagulated.
3. Connect plug(s) to appropriate bipolar electrosurgical generator.
4. Immerse distal tip of probe into saline, then remove. Immediately after removing tip from water source, activate generator to test probe.  
**Note:** Water should bubble on tip of probe. If an abnormality is detected which would prohibit proper working condition, do not use. Turn generator off.
5. Immediately prior to insertion of device into endoscope, immerse tip of probe in saline or distilled water for at least 5 seconds to activate Hydromer hydrophilic grafting.
6. Introduce device into endoscope and advance in small increments to desired site.
7. Following electrosurgical unit manufacturer's instructions for settings, verify desired settings. Verify proper position of probe, then activate electrosurgical unit. **Warning:** Contact of probe with endoscope during electrosurgery may cause grounding which could result in injury to patient and/or operator as well as damage to endoscope and/or device.
8. Proceed with bipolar coagulation.
9. Upon completion of procedure, turn electrosurgical unit off. Remove device from endoscope and disconnect from electrosurgical unit. Dispose of probe per institutional guidelines for biohazardous medical waste.

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Hydromer is a registered trademark of Hydromer Inc.
Quicksilver is a trademark of Cook Endoscopy.
Cook is a registered trademark of Cook Incorporated.
**Cook Ireland Ltd.**