**Memory™ Hard Wire Baskets**

**Cestos Rígidos Memory™**

**Cestas Memory™ Alambre Duro**

18939/0205

**COOK**® Wilson-Cook Medical  
GI Endoscopy

English . . . . .	1
Português/Portuguese . . . . .	5
Español/Spanish . . . . .	10

## **INTENDED USE**

The Memory Hard Wire Basket is used for the endoscopic removal of biliary stones and foreign bodies. This device is supplied sterile and is reusable if the device integrity is intact.

## **NOTES**

Do not use this device for any purpose other than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device.

If the product package is open or damaged when received, do not use this device.

Wilson-Cook devices should be stored in a dry location, away from temperature extremes.

## **CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic procedure to be performed in conjunction with the removal of the foreign body.

Perforation.

If the Memory Hard Wire Basket is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

## **POTENTIAL COMPLICATIONS**

Potential complications associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Other potential complications associated with basket extraction include, but are not limited to: impaction of the object, aspiration of foreign body, localized inflammation, pressure necrosis.

## **PRECAUTIONS**

The coordination of endoscope accessory channel size with compatible devices is essential in order to obtain optimal results during a procedure. The sheath of the Memory Hard Wire Basket is 7.0 FR and requires a minimum channel size of 2.8 mm.

If this device is to be used for the removal of biliary stones, assessment of the stone size and ampullary orifice must be made to determine the necessity of sphincterotomy. In the event sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

Injection of contrast during ERCP should be monitored fluoroscopically.

Overfilling of the pancreatic duct may cause pancreatitis.

Memory Hard Wire Baskets should never be coiled in less than an 8-inch (20 cm) diameter.

## **WARNINGS**

Wilson-Cook Memory Hard Wire Baskets are **not** compatible with the Soehendra® Lithotripter or any other mechanical lithotripter.

When retrieving a foreign body from the upper GI tract, extreme caution must be exercised to avoid occlusion of the airway.

Surgical intervention may be required if impaction occurs.

## **INSTRUCTIONS FOR USE**

1. Upon removal of the device from the package and prior to each use, uncoil the device, extend the basket from the sheath and visually inspect with particular attention to kinks, bends, or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Wilson-Cook for return authorization if an abnormality is noted upon initial inspection of the device.
2. With the basket fully retracted into the sheath, insert the device into the endoscope accessory channel. Advance the device through the channel, in short increments, until the basket sheath exits the endoscope.
3. Endoscopically or fluoroscopically visualize the extraction target and advance the basket to the appropriate position. Fluoroscopic visualization may be enhanced by injecting contrast via the side arm fitting. Place the sheath of the basket beyond the target before extending the basket to entrap the extraction target.
4. After confirming the desired position of the basket sheath relative to the target, advance the basket out of the sheath by pushing forward on the handle. Place the basket around the extraction target and slowly retract the handle until entrapment is achieved. **Note:** It may be necessary to move the basket back and forth along the target in order to capture the object.
5. Withdraw the basket into the channel while maintaining entrapment. Support the foreign body or biliary stone against the endoscope tip, and withdraw the endoscope from the patient. Maintain endoscopic visualization to ensure continued entrapment of the object. While withdrawing the device from the endoscope, wipe the secretions from the outer sheath. **Caution:** If difficulty is encountered when removing the basket from the duct, moderate force may be applied by pulling on the handle to manually fracture the stone. If passage is still restricted, surgical intervention may be necessary.
6. Prepare or dispose of the extracted specimen per institutional guidelines.

7. Upon completion of the procedure, clean and prepare the device for future use following the cleaning instructions below. During cleaning, inspect the integrity and function of the device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use. Dispose of the device per institutional guidelines for biohazardous medical waste.

## **CLEANING INSTRUCTIONS**

**Caution:** During handling, coil the device in a minimum 8-inch diameter to prevent product damage.

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to the Luer lock fitting on the device, then irrigate.
2. Extend the basket, then submerge the solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
3. Beginning at the handle and working toward the tip, wipe the device with a textured cloth soaked in enzymatic solution.
4. Irrigate the inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to the Luer lock fitting on the handle of the device. Repeat this step until all visible secretions are removed.
5. Gently brush the basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
6. Remove the device from the enzymatic solution, then thoroughly rinse the inner sheath using a syringe filled with clean water. Rinse the handle, outer sheath and extended basket thoroughly with clean, running water.
7. Ultrasonically clean the device, with the basket extended, for a minimum of 15 minutes at 25°C (77°F) or per the ultrasonic cleaner manufacturer's instructions.
8. Remove the device from the ultrasonic cleaner, then thoroughly rinse the inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse the external portions of the device and the basket with clean, running water.
9. Gently wipe the external portions of the device dry with a lint-free cloth. Air dry the lumen and the basket by hanging the device, with the basket extended, for three hours, or by using medical grade (oil free) forced air through the Luer lock fitting on the handle.
10. Package the device for sterilization following AAMI recommended practices. **Note:** AAMI stands for Association for the Advancement of Medical Instrumentation.
11. Sterilize the device using the set points for the desired sterilization method shown in the tables that follow.

## ETO CYCLE PARAMETERS

### PARAMETER

Prevacuum Pressure  
Relative Humidity  
Humidity Dwell  
Chamber Temperature  
Exposure Pressure  
Exposure Dwell  
ETO Concentration  
Post Vacuums  
Post Vacuum Pressure  
Air Wash<sup>2</sup>  
Aeration Temperature  
Aeration Dwell  
Air Changes During Aeration

### SET POINT

2.0 psia (-25.9" Hg Vac)  
60% RH  
20 minutes  
54.4°C (130°F)  
To Be Determined<sup>1</sup>  
2 hours  
600 mg/L  
3  
2.0 psia (-25.9" Hg Vac)  
5 minutes  
54.4°C (130°F)  
12 hours  
90 per hour

<sup>1</sup> To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

<sup>2</sup> Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

**Note:** ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilizing, it should be taken out of service immediately.

## AUTOCLAVE CYCLE PARAMETERS

### PARAMETER

Cycle Type  
Prevacuum Pressure  
Chamber Exposure Temperature  
Chamber Exposure Pressure  
Exposure Dwell

### SET POINT

Prevacuum  
2.5 psia (-24.8" Hg Vac)  
132.2°C (270°F)  
41.8 psia (27.1 psig)  
5 minutes

**Note:** Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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Memory is a trademark of Wilson-Cook Medical Inc.
Soehendra is a registered trademark of Wilson-Cook Medical Inc.
Cook is a registered trademark of Cook Incorporated.
**Wilson-Cook Medical Inc.**
**Cook Ireland Ltd.**