COOK®

MEDICAL

PreCurved Tip ERCP Catheters

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INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

NOTES

This device is reusable if integrity is intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Reusability of a device depends, in large part, on care of device by user. Factors involved in prolonging life of this device include, but are not limited to: advancing through endoscope accessory channel in short increments, gently withdrawing device from channel, avoiding loops and/or kinks in catheter, coiling device in a minimum 8-inch diameter, thorough cleaning following instructions included in this booklet.

INSTRUCTIONS FOR USE

1. Gently advance stylet completely into catheter. Lock stylet cap into catheter fitting.
2. Introduce catheter into endoscope and advance in small increments until device is endoscopically visualized exiting scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm position of device by attaching a syringe to side port of catheter fitting.
4. To insert a compatible wire guide, remove stylet from catheter and advance wire guide prior to removing catheter from endoscope. **Note:** Before each introduction of a wire guide into a device, wet wire guide in a sterile water or saline bath. A contrast-filled catheter may make wire guide difficult to advance. Flush contrast from catheter with sterile water or saline before each introduction of wire guide, observing necessary cautions if device is positioned in pancreatic duct.

5. Upon completion of procedure, clean and prepare catheter for future use following cleaning instructions.

CLEANING INSTRUCTIONS

1. While withdrawing catheter from endoscope, wipe excess secretions from outer sheath.
2. Immediately after use, remove stylet wire from catheter and attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting. Block open port on catheter fitting, then irrigate inner lumen of catheter.
3. Submerge solution-filled catheter and stylet in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
4. Beginning at Luer lock fitting and working toward tip, wipe catheter with a textured cloth soaked in enzymatic solution. Beginning at hub and working downward, wipe stylet.
5. Block open port on catheter fitting, then irrigate inner lumen of catheter using a syringe containing a minimum of 20 cc of enzymatic solution attached to Luer lock fitting. Repeat this step until all visible secretions are removed.
6. Remove catheter and stylet from enzymatic solution, then thoroughly rinse inner lumen of catheter using a syringe filled with clean water. Rinse stylet and outer portions of catheter thoroughly with clean, running water.
7. Ultrasonically clean catheter and stylet for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove catheter and stylet from ultrasonic cleaner. To remove all residue, thoroughly rinse inner lumen of catheter using a syringe containing a minimum of 20 cc of clean water. Thoroughly rinse stylet and outer portions of catheter with clean, running water.
9. Gently wipe external portions of stylet and catheter dry with a lint-free cloth. Hang device to air dry inner lumen of catheter, or use medical grade (oil free) forced air through Luer lock fitting.
10. Insert stylet into catheter, leaving approximately 1 cm extending from Luer lock fitting.
11. Package device for sterilization following AAMI recommended practices.
12. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER

Prevacuum Pressure
Relative Humidity
Humidity Dwell
Chamber Temperature
Exposure Pressure
Exposure Dwell
ETO Concentration
Post Vacuums
Post Vacuum Pressure
Air Wash²
Aeration Temperature
Aeration Dwell
Air Changes During Aeration

SET POINT

2.0 psia (25.9" Hg Vac)
60% RH
20 minutes
54.4°C (130°F)
To Be Determined¹
2 hours
600 mg/L
3
2.0 psia (25.9" Hg Vac)
5 minutes
54.4°C (130°F)
12 hours
90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilizing, it should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER

Cycle Type
Prevacuum Pressure
Chamber Exposure Temperature
Chamber Exposure Pressure
Exposure Dwell

SET POINT

Prevacuum
2.5 psia (24.8" Hg Vac)
132.2°C (270°F)
41.8 psia (27.1 psig)
5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No. 12 - 1994.

"Sterilization, Part 1 - Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

Cook is a registered trademark of Cook Incorporated.

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

RxOnly

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland

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