Cunningham-Cotton Sleeve

Cunningham-Cotton 套管

Cunningham-Cotton manchet

Cunningham-Cotton huls

Cathéter Cunningham-Cotton

Cunningham-Cotton-Hülse

Χιτώνιο Cunningham-Cotton

Catetere Cunningham-Cotton

Cateter Cunningham-Cotton

Manguito Cunningham-Cotton

Cunningham-Cotton hölje

18966/0305

COOK®

Endoscopy

English.	1
Chinese/中文	2
Danish/Dansk	3
Dutch/Nederlands	5
French/Français	7
German/Deutsch	9
Greek/Ελληνικά	11
Italian/Italiano	13
Portuguese/Português	15
Spanish/Español	17
Swedish/Svenska	19

INTENDED USE

This device is used for endoscopic cannulation of the ductal system. The outer sleeve allows maintained access to the biliary duct to facilitate papillotomy/sphincterotomy, biliary stent placement and/or nasal biliary drainage. This device is supplied sterile and intended for single use only.

NOTE

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Endoscopy for return authorization.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

Cook Endoscopy wire guides are recommended for use with this device.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

1. Introduce and advance the device over a prepositioned wire guide.
2. Under fluoroscopic monitoring, advance the device to the desired location. **Note:** This device is equipped with two radiopaque bands on the dilator portion of the device and one radiopaque band on the sleeve portion of the device in order to fluoroscopically visualize the position of the catheter in the duct.
3. Remove the dilator portion of the device and the wire guide (if necessary), then fluoroscopically confirm proper position of the device.
4. Proceed with the exchange of other compatible Cook Endoscopy devices through the sleeve portion of the device.
5. Upon completion of the procedure, remove the device from the accessory channel of the endoscope and dispose of per institutional guidelines for biohazardous medical waste.

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© 2005 Cook Endoscopy
Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.