Esophageal Dilation Balloon

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Dilatationsballon til øsofagus

Oesofageale dilatatieballon

Ballonnet de dilatation oesophagienne

Ösofagusdilatationsballon

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INTENDED USE

This device is used to endoscopically dilate strictures of the esophagus.

PRODUCT DESCRIPTION

This balloon can be inflated to three distinct diameters, as indicated on package and catheter tag. This device is used in conjunction with an inflation device.

NOTES

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Endoscopy for return authorization.

Do not use for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications to dilation include, but are not limited to: asymptomatic rings, webs, or strictures; known or suspected perforation; severe inflammation or scarring near the dilation site.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for the minimum channel size requirement for this device.

Do not pre-inflate the balloon.

WARNINGS

During dilation, do not inflate balloon beyond the maximum indicated inflation pressure.

The balloon must be inflated with water only.

Do not advance the balloon dilator if resistance is encountered. Assess the cause of resistance to determine if dilation should be reattempted.

INSTRUCTIONS FOR USE

1. Attach the balloon to a 60 ml (cc) inflation device with gauge to monitor balloon pressure.
2. To facilitate passage through the endoscope, apply negative pressure to the catheter.
3. Remove the protective sheath from the balloon.
4. Apply a lubricating agent to the balloon to facilitate passage through the endoscope accessory channel.

5. Maintain balloon deflation with negative pressure and introduce into the accessory channel of the endoscope, advancing in short increments until the dilator is completely visualized endoscopically. Monitor endoscopically until the balloon is in the desired position within the stricture.
6. Inflate the balloon to the pressure corresponding to the smallest balloon diameter and maintain until desired dilation is achieved. To achieve increasingly larger balloon diameters, increase pressure as indicated on catheter tag. **Note:** Do not exceed the maximum indicated inflation pressure.
7. Apply a vacuum and remove all fluid from the balloon while observing the balloon endoscopically.
8. Remove the deflated balloon from the accessory channel.
Caution: A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.
9. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.