COOK®

MEDICAL

Biliary Wire Guide (ES)

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INTENDED USE

This device is used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: perforation, allergic reaction to contrast or medication.


PRECAUTIONS

Use of this wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

This wire guide may be left in place when used under normal conditions in conjunction with, and according to the instruction for, a compatible Cook electro-surgical (ES) device.

INSTRUCTIONS FOR USE

1. Prior to removing wire guide from holder, flush with 30 cc of sterile water. (See fig.)
2. Flush endoscope accessory channel and/or lumen of device with sterile water, **then insert wire guide floppy end first. Note:** For best results, wire guide should be kept wet, if applicable.
3. Fluoroscopically monitor wire guide advancement in ductal system.
4. Proceed with exchange of compatible wire-guided accessories over wire guide.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.