Surgisis[®] Biodesign[™]
URETHRAL SLING

FP0007-3B

Manufacturer:

Cook Biotech Incorporated
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.
Phone: (812) 339-2235
Toll Free: (800) 457-4500
www.cookmedical.com

EU Representative:

Cook Ireland, Ltd.
O' Halloran Road
National Technological Park
Limerick, Ireland
Tel: +353 61 334440
Fax: +353 61 334441

© COPYRIGHT COOK BIOTECH 2007

ENGLISH

INTENDED USE:

The Cook[®] Surgisis[®] Biodesign[™] Urethral Sling is used to provide pubourethral support for the treatment of urinary incontinence resulting from hypermobility or intrinsic sphincter deficiency. The Urethral Sling is intended for implantation to reinforce soft tissues where weakness exists.

This device is supplied sterile in peel-open packages and is intended for one-time use.

[Rx ONLY] This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS: This device is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

PRECAUTIONS:

- **Do not resterilize.** Discard all open and unused portions of the Urethral Sling.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that device is rehydrated prior to securing the device to the surgical site with sutures, staples, or other fasteners.
- Device performance has not been evaluated with suture spacing greater than 2 mm.
- Ensure that all layers of the Urethral Sling are secured when fastening the device to the surgical site.
- Excess tension on the device may cause urinary retention in the device recipient when being treated for urinary incontinence.

POTENTIAL COMPLICATIONS: The following complications are possible with the use of surgical graft materials. If any of these conditions occur, the device should be removed.

- Infection
- Acute or Chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

STORAGE: This device should be stored in a clean, dry location at room temperature.

STERILIZATION: This device has been sterilized with ethylene oxide. Do not resterilize.

SUGGESTED INSTRUCTIONS FOR USING THE URETHRAL SLING

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Always handle sheets using aseptic technique, minimizing contact with latex gloves.

Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 50 ml of room temperature sterile saline or sterile lactated Ringer's solution for each Urethral Sling

1. Using aseptic technique, remove the the Urethral Sling inner pouch from its outer bag, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the Urethral Sling with the sterile forceps.
3. Place the Urethral Sling into the sterile dish in the sterile field.
4. Add to the dish at least 50 ml of the rehydration fluid for each Urethral Sling.
5. Allow the Urethral Sling to rehydrate for at least ten minutes.
6. Prepare the surgical site using standard surgical techniques.
7. Using aseptic technique, trim the Urethral Sling to fit the surgical site. (Note: An alternative method is to cut the Urethral Sling to size prior to rehydration. If this method is selected, be sure to rehydrate the Urethral Sling prior to securing the device into place. See step 5.)
8. Using aseptic technique, transfer the Urethral Sling to the surgical site and secure the device into place, avoiding excess tension. (Note: Surgical experience indicates that securing the Urethral Sling with close tissue approximation produces better outcomes.)
9. Complete the surgical procedure.
10. Discard any unused portions of the Urethral Sling according to institution guidelines for medical waste.

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