**Biodesign®
Surgisis™ 8-Layer Tissue Graft
Surgisis™ Hernia Graft**

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CE 0123

INTENDED USE

The Cook® Biodesign™ Surgisis® Hernia Graft and the 8-Layer Tissue Graft are intended to be implanted to reinforce soft tissues where weakness exists. Indications for use include the repair of a hernia or body wall defect. Biodesign Surgisis is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This device is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

PRECAUTIONS

- **Do not resterilize.** Discard all open and unused portions of the device.
- The Biodesign Surgisis is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that device is rehydrated prior to cutting, suturing, stapling, tacking or loading of the device laparoscopically.
- Ensure that all layers of the Biodesign Surgisis are secured when suturing or stapling.
- Place device in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing or stapling more than one device together may decrease device performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the Biodesign Surgisis.
- Extended rehydration or excessive handling could lead to partial delamination of superficial layers of the Biodesign Surgisis.
- Biodesign Surgisis has been reported to be safe in clean and clean-contaminated hernia repair.¹ However, care should be exercised when it is placed in critically ill patients or in severely contaminated wounds.

USE OF ANTIMICROBIALS

Because the Biodesign Surgisis is at times used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.² In these cases both antibiotic prophylaxis of the patient and antimicrobial soaking of the graft have been used. Typical flora can be expected to include a variety of aerobic and facultative anaerobic organisms, including, but not limited to, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.
- Antibacterial prophylaxis, if chosen, should be started prior to surgery and continued post-operatively.²

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the Biodesign Surgisis. For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,^{3,2} while povidone iodine,⁴ bacitracin,¹³ polymyxin B,⁷ and vancomycin⁸ have all been reported to slow or inhibit wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with the Biodesign Surgisis with respect to hernia repair.

POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of any prosthesis may include, but are not limited to, infection, inflammation, adhesion, fistula formation, seroma formation, hematoma, and recurrence of tissue defect.

Complications, such as delayed wound infection, hernia recurrence, and the need for re-operation, should be reasonably expected in patients who are critically ill or who have severely contaminated abdomens.

STORAGE

This device should be stored in a clean, dry location at room temperature.

STERILIZATION

This device has been sterilized with ethylene oxide.

SUGGESTED INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Always handle the Biodesign Surgisis using aseptic technique, minimizing contact with latex gloves.

Required Materials

- Sterile forceps
- Rehydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.

1. Using aseptic technique, remove the sterile inner package from the outer packaging and place it in the sterile field.

- Open the sterile inner package carefully. Aseptically, remove the Biodesign Surgisis from the inner package.
- Rehydrate the device in room temperature, sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. For the Hernia Graft, a rehydration time of greater than 1 minute is not required. For the 8-Layer Tissue Graft, rehydrate for at least five (5) minutes prior to cutting, suturing, stapling, tacking or loading the device laparoscopically. If the device is provided in a tray, it can be rehydrated in the tray. (See Use of Antimicrobials)
- Prepare the hernia repair site using standard surgical techniques.
- Using aseptic technique, trim the Biodesign Surgisis to fit the site, providing an allowance for overlap. Fundamental surgical principles suggest recurrence can be minimized if the mesh overlaps the surrounding tissue by at least 4-5 cm in all directions.⁹

NOTE: If the device is cut too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original tissue defect or development of a defect in the adjacent tissues.

- Using aseptic technique, transfer the Biodesign Surgisis to the surgical site and suture or staple into place, avoiding excess tension.

NOTE: Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

NOTE: Surgical experience indicates that suturing or stapling the Biodesign Surgisis with close tissue approximation produces better outcomes. Fundamental surgical principles suggest a suture spacing approximately equal to suture bite depth.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

NOTE: The Hernia Graft has been designed with an open matrix structure to optimize tissue incorporation and angiogenesis.

- Complete the standard surgical procedure.
- Discard any unused portions of the Biodesign Surgisis following standard technique for disposal of medical waste.

REFERENCES

- Helton WS, et al., *Short-term Outcomes with Small Intestinal Submucosa for Ventral Abdominal Hernia* Arch Surg. 2005; 140: 549-562
- Mangram AJ, et al., *Guideline for the prevention of surgical site infection, 1999.* Centers for Disease Control and Prevention, 1999.
- Bang K, et al., *Scand J Plast Reconstr Surg Hand Surg* 1998;32:129-33.
- Petroutsos G, et al., *Arch Ophthalmol* 1983;101:1775-8.
- Nelson JD, et al., *Curr Eye Res* 1990;9:277-85.
- Kjølseth D, et al., *J Am Coll Surg* 1994;179:305-12.
- Nakamura M, et al., *Curr Eye Res* 1993;12:733-40.
- Petroutsos G, et al., *Int Ophthalmol* 1984;7:65-9.
- Klinge U, et al., *World J Surg* 2005;29:1066-72.

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LITTERATUR

- 1 Helton WS, et al., *Short-term Outcomes with Small Intestinal Submucosa for Ventral Abdominal Hernia* Arch Surg. 2005; 140: 549-562
- 2 Mangram AJ, et al., *Guideline for the prevention of surgical site infection*, 1999. Centers for Disease Control and Prevention, 1999.
- 3 Bang K, et al., *Scand J Plast Reconstr Surg Hand Surg* 1998;32:129-33.
- 4 Petroustos G, et al., *Arch Ophthalmol* 1983;101:1775-8.
- 5 Nelson JD, et al., *Curr Eye Res* 1990;9:277-85.
- 6 Kjolseth D, et al., *J Am Coll Surg* 1994;179:305-12.
- 7 Nakamura M, et al., *Curr Eye Res* 1993;12:733-40.
- 8 Petroustos G, et al., *Int Ophthalmol* 1984;7:65-9.
- 9 Klinge U, et al., *World J Surg* 2005;29:1066-72.

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QUELLEN

- 1 Helton WS, et al., *Short-term Outcomes with Small Intestinal Submucosa for Ventral Abdominal Hernia* Arch Surg. 2005; 140: 549-562
- 2 Mangram AJ, et al., *Guideline for the prevention of surgical site infection, 1999.* Centers for Disease Control and Prevention, 1999.
- 3 Bang K, et al., *Scand J Plast Reconstr Surg Hand Surg* 1998;32:129-33.
- 4 Petroutsos G, et al., *Arch Ophthalmol* 1983;101:1775-8.
- 5 Nelson JD, et al., *Curr Eye Res* 1990;9:277-85.
- 6 Kjolseth D, et al., *J Am Coll Surg* 1994;179:305-12.
- 7 Nakamura M, et al., *Curr Eye Res* 1993;12:733-40.
- 8 Petroutsos G, et al., *Int Ophthalmol* 1984;7:65-9.
- 9 Klinge U, et al., *World J Surg* 2005;29:1066-72.

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- 1 Helton WS, et al., *Short-term Outcomes with Small Intestinal Submucosa for Ventral Abdominal Hernia* Arch Surg. 2005; 140: 549-562
- 2 Mangram AJ, et al., *Guideline for the prevention of surgical site infection*, 1999. Centers for Disease Control and Prevention, 1999.
- 3 Bang K, et al., *Scand J Plast Reconstr Surg Hand Surg* 1998;32:129-33.
- 4 Petroutsos G, et al., *Arch Ophthalmol* 1983;101:1775-8.
- 5 Nelson JD, et al., *Curr Eye Res* 1990;9:277-85.
- 6 Kjolseth D, et al., *J Am Coll Surg* 1994;179:305-12.
- 7 Nakamura M, et al., *Curr Eye Res* 1993;12:733-40.
- 8 Petroutsos G, et al., *Int Ophthalmol* 1984;7:65-9.
- 9 Klinge U, et al., *World J Surg* 2005;29:1066-72.

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REFERENSER

- 1 Helton WS, et al., *Short-term Outcomes with Small Intestinal Submucosa for Ventral Abdominal Hernia* Arch Surg. 2005; 140: 549-562
- 2 Mangram AJ, et al., *Guideline for the prevention of surgical site infection, 1999.* Centers for Disease Control and Prevention, 1999.
- 3 Bang K, et al., *Scand J Plast Reconstr Surg Hand Surg* 1998;32:129-33.
- 4 Petroustos G, et al., *Arch Ophthalmol* 1983;101:1775-8.
- 5 Nelson JD, et al., *Curr Eye Res* 1990;9:277-85.
- 6 Kjolseth D, et al., *J Am Coll Surg* 1994;179:305-12.
- 7 Nakamura M, et al., *Curr Eye Res* 1993;12:733-40.
- 8 Petroustos G, et al., *Int Ophthalmol* 1984;7:65-9.
- 9 Klinge U, et al., *World J Surg* 2005;29:1066-72.

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.
Phone: 812 339-2235
Toll Free: 800 457-4500

EU REPRESENTATIVE

COOK IRELAND
O'Halloran Road
National Technological Park
Limerick, IRELAND
Phone: +353 61 334440

COOK (CANADA) INC.
111 A Drive
Stouffville, Ontario
L4A 7X5 CANADA
Phone: 905 640-7110
Toll Free: 800 668-0300

WILLIAM A. COOK
AUSTRALIA PTY. LTD.
Brisbane Technology Park
12 Electronics Street
Eight Mile Plains
Brisbane, QLD 4113 Australia
Phone: +61 7 38 41 11 88

www.cookmedical.com

© COPYRIGHT COOK BIOTECH 2008 SEPTEMBER/2008