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U _ C L S T _ R E V . 1

SUGGESTED CLEANING AND STERILIZATION INSTRUCTIONS

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

There are no known contraindications.

CAUTION: This product must be cleaned immediately after use. Blood, mucus and other debris, if allowed to dry, can render the device inoperable.

1. Remove all residual organic matter, blood and irrigation solution to allow the sterilization or disinfection medium to contact the surface of the device. **NOTE: Ethylene oxide (ETO) gas and disinfection solutions do not effectively penetrate residual organic matter.**
2. Clean all surfaces with a solution of mild detergent and water. Use a soft brush to remove all residual organic matter, blood and irrigation solution from all surfaces.
3. Ultrasonic cleaning of the components may be performed. Ultrasonic cleaning should only be performed with non-abrasive and non-corrosive solutions.
4. Rinse thoroughly with clean water.
5. Dry all exterior surfaces with a soft cloth or an air pistol.
6. Dry all openings with an air pistol.

NOTE: Prior to sterilization and reuse, inspect device to verify that it functions properly. All edges and surfaces should be smooth and free of abrasions.

ETHYLENE OXIDE GAS STERILIZATION

Routine ethylene oxide (ETO) gas sterilization is the sterilization method of choice.

ETO gas sterilization and aeration should be performed in accordance with AAMI Standards and Recommended Practices with a validated sterilization cycle.

AUTOClave STERILIZATION

This device can be sterilized by steam autoclave. Autoclave sterilization should be performed in accordance with AAMI Good Hospital Practices with a validated sterilization cycle. Autoclave temperature should not exceed 275°F.

DISINFECTION

Chemical sterilization should be performed in accordance with AAMI Good Hospital Practices.

NOTE: Only EPA registered hard surface disinfectants capable of high-level disinfection should be used.

1. Carefully follow the disinfectant manufacturer instructions.
2. Devices to be disinfected should be placed in a plastic basin to avoid scratches and electrolytic corrosion.
3. Rinse disinfected devices thoroughly to remove all traces of residue.

AAMI/ISO 11135-1994-Medical Devices-Validation and Routine Control of Ethylene Oxide Sterilization.

ANSI/AAMI ST79:2006-Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

ANSI/AAMI ST41:1999-Ethylene Oxide Sterilization in Health Care Facilities: Safety and Effectiveness.

ANSI/AAMI ST58:2005-Chemical Sterilization and High-Level Disinfection in Health Care Facilities.

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ANSI/AAMI ST79:2006-Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities (卫生保健机构蒸汽灭菌及灭菌保证综合指南)。

ANSI/AAMI ST41:1999-Ethylene Oxide Sterilization in Health Care Facilities: Safety and Effectiveness (卫生保健机构环氧乙烷灭菌：安全及有效性)。

ANSI/AAMI ST58:2005-Chemical Sterilization and High-Level Disinfection in Health Care Facilities (卫生保健机构化学灭菌及强力消毒说明)。

COOK UROLOGICAL, INC.
COOK OB/GYN
1100 West Morgan Street
Spencer, Indiana 47460, USA

COOK IRELAND LTD.
O'Halloran Road
National Technology Park
Limerick, Ireland

© COPYRIGHT COOK 2008

200804