**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

**Aspiracath™**

U \_ J - A S P \_ R E V . 2

EN  
3  
CS  
3  
DA  
4  
DE  
4  
EL  
5  
ES  
5  
FR  
6  
HU  
6  
IT  
7  
NL  
7  
NO  
8  
PL  
8  
PT  
9  
SV  
10  
ZH  
10

## ASPIRACATH™

Used to aspirate cervical mucous. Intended for one-time use.

### DESCRIPTION

- Polyethylene aspiration catheter with integrated nylon plunger

**CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.**

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

### SUGGESTED INSTRUCTIONS FOR USING ASPIRACATH™

#### CONTRAINDICATIONS

Use of this device is contraindicated in the presence of pregnancy, cervical cancer, purulent infections in the vagina or cervix, untreated uterine anomaly.

1. Place the patient in lithotomy position.
2. With the plunger fully inserted in the catheter, insert the device into the cervix which is approximately 1-3 cm.

**NOTE: Do not use undue force while advancing the catheter if resistance is met.**

3. Gently withdraw the plunger from the catheter while removing the device from the cervix.
4. Advance the plunger into the catheter to expel the mucous contents onto a glass slide. Allow sample to air dry and it is then examined microscopically.
5. Dispose of device.

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**COOK IRELAND LTD.**