Embryo Transfer Catheter

Katetr pro embryotransfer

Embryo transfereringskateter

Embryotransferkatheter

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U _ K - J - E T C _ R E V . 1

EMBRYO TRANSFER CATHETER

Used to place *in vitro* fertilized (IVF) embryos into the uterine cavity. Intended for one-time use.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

This catheter should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent pregnancy (or is currently pregnant), or if the patient currently has an intrauterine device.

PRECAUTIONS

The assisted reproduction catheter (and any other accessories used during this procedure) should be comprised of embryo compatible materials.

Infection may occur due to bacterial contamination of the device during vaginal manipulation, and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing the catheter (and any other accessories used) with sterile, compatible culture media, and closely adhering to sterile techniques.

Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.

SUGGESTED INSTRUCTIONS FOR USING EMBRYO TRANSFER CATHETER

1. Place patient in a lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. Use wire guide or malleable/flexible obturator with introducing/guide catheter/cannula to obtain placement if necessary. Remove wire guide or malleable/flexible obturator prior to placement of inner/delivery catheter in introducing/guide catheter or for the locking version, lock the inner/delivery catheter into the introducing/guide catheter prior to passing into the patient.
4. A syringe is attached to the inner/delivery catheter and embryos are aspirated into the end of the catheter.
5. The inner/delivery catheter, previously loaded with embryos, is then introduced directly into the uterine cavity or through the introducing/guide catheter into the uterine cavity.
6. Embryos are then expelled into the uterine cavity.
7. Remove the catheter/catheters, check for retention of embryos, and discard.

NOTE: One cell mouse embryo tested and passed with 75% or greater blastocyst rate. USP endotoxin (LAL) tested and passed with 20 EU's or less per device. Testing is conducted on a lot-to-lot (batch) basis.

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COOK IRELAND LTD.