Intratubal Transfer Set

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U _ K - J - I T S _ R E V . 1

INTRATUBAL TRANSFER SET

Used to inject either sperm, gametes, or embryos into the uterine ostium of the fallopian tube via ultrasound guidance. Intended for one-time use.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

This catheter set should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent pregnancy (or is currently pregnant), or if the patient currently has an intrauterine device.

PRECAUTIONS

The assisted reproduction catheter (and any other accessories used during this procedure) should be comprised of embryo compatible materials.

Infection may occur due to bacterial contamination of the catheter during vaginal and transcervical manipulation, and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing the catheter (and any other accessories used) with sterile, compatible culture media, administering prophylactic antibiotics, and closely adhering to sterile techniques.

Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.

SUGGESTED INSTRUCTIONS FOR USING INTRATUBAL TRANSFER SET

NOTE: It is important to determine the exact catheter volume prior to use. To determine catheter volume, aspirate transfer media or other appropriate liquid into the catheter, utilizing a calibrated syringe.

1. Under ultrasound guidance, pass the preassembled malleable obturator and guiding catheter transvaginally through the uterus and into the tubal ostium.

NOTE: The obturator can be withdrawn to manipulate the preset curve on the catheter and facilitate placement.

2. Once in position, remove the malleable obturator.
3. If desired, the wire guide may be introduced into the guiding catheter to facilitate cannulation of the fallopian tube.
4. Introduce the preloaded delivery catheter into the guiding catheter, advance the delivery catheter approximately 2 cm into the fallopian tube and instill the desired samples.
5. Remove the device and discard.

NOTE: One cell mouse embryo tested and passed with 75% or greater blastocyst rate. USP endotoxin (LAL) tested and passed with 20 EU's or less per device. Testing is conducted on a lot-to-lot (batch) basis.

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COOK IRELAND LTD.