T.V.P. Embryo Transfer Catheter

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U _ T V P C _ R E V . 1

T.V.P. EMBRYO TRANSFER CATHETER

Used to place *in vitro* fertilized (IVF) embryos into the uterine cavity. Intended for one-time use.

CONTRAINDICATIONS

This catheter should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent pregnancy (or is currently pregnant), or if the patient currently has an intrauterine device.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

PRECAUTIONS

The assisted reproduction catheter (and any other accessories used during this procedure) should be comprised of embryo compatible materials.

Infection may occur due to bacterial contamination of the device during vaginal manipulation, and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing the catheter (and any other accessories used) with sterile, compatible culture media, and closely adhering to sterile techniques.

Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.

SUGGESTED INSTRUCTIONS FOR USING EMBRYO TRANSFER CATHETER WITH TRANSVAGINAL PROBE

1. Place the patient in lithotomy position.
2. Perform a vaginal ultrasound and measure the endocervical length, the endometrial length, and the cervico-uterine angle.
3. Insert open sided vaginal speculum to expose the cervix.
4. Wash the vagina with saline, aspirate the cervical mucus and wash the endocervix with media.
5. Load embryos into the inner delivery catheter. Lock the delivery catheter into the introducing/guide catheter.
6. The stopper on the catheter is placed at a distance of one centimeter plus the length of the endocervix, and an appropriate curvature is placed on the catheter consistent with the previously measured cervico-uterine angle.
7. The catheter is then introduced into the uterine cavity, the open-sided speculum is removed, the vaginal ultrasound probe is introduced into the vagina. The catheter is visualized, the lock is released and the inner catheter is advanced to within 1-2 cm of the uterine fundus.
8. In difficult cases use wire guide or malleable/flexible obturator with introducing/guide catheter/cannula to obtain placement if necessary. Remove wire guide or malleable/flexible obturator prior to placement of inner/delivery catheter.
9. Embryos are then expelled into the uterine cavity.
10. Remove the catheter/catheters, check for retention of embryos, and discard.

NOTE: One cell mouse embryo tested and passed with 75% or greater blastocyst rate. USP Endotoxin (LAL) tested and passed with 20 EU's or less per device. Testing is conducted on a lot-to-lot (batch) basis.

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COOK IRELAND LTD.