EN
3

Bakri Postpartum Balloon

Instructions for Use

CS
5

Tampónovací balónek Bakri

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DA
8

Bakri postpartum-ballon

Brugsanvisning

DE
11

Bakri Postpartum-Ballon

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EL
13

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

BAKRI POSTPARTUM BALLOON

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

This device is intended to provide temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.

CONTRAINDICATIONS

- Arterial bleeding requiring surgical exploration or angiographic embolization
- Cases indicating hysterectomy
- Pregnancy
- Cervical cancer
- Purulent infections in the vagina, cervix, or uterus
- Untreated uterine anomaly
- Disseminated intravascular coagulation
- A surgical site that would prohibit the device from effectively controlling bleeding

WARNINGS

- This device is intended as a temporary means of establishing hemostasis in cases indicating conservative management of postpartum uterine bleeding.
- **The Bakri Postpartum Balloon is indicated for use in the event of primary postpartum hemorrhage within 24 hours of delivery.**
- The device should not be left indwelling for more than 24 hours.
- **The balloon should be inflated with a sterile liquid such as sterile water, sterile saline, or lactated ringers solution. The balloon should never be inflated with air, carbon dioxide or any other gas.**
- The maximum inflation is 500 ml. Do not overinflate the balloon. **Overinflation of the balloon may result in the balloon being displaced into the vagina.**
- Clinical data to support the safety and effectiveness of the Bakri Postpartum Balloon in the setting of uterine atony are limited. Patients in whom this device is being used should be closely monitored for signs of worsening bleeding and/or disseminated intravascular coagulation (DIC). In such cases, emergency intervention per hospital protocol should be followed.
- There are no clinical data to support use of this device in the setting of DIC.
- Patient monitoring is an integral part of managing postpartum hemorrhage. Signs of deteriorating or non-improving condition should lead to a more aggressive treatment and management of patient uterine bleeding.
- Patient urine output should be monitored while the Bakri Postpartum Balloon is in use.

PRECAUTIONS

Avoid excessive force when inserting the balloon into the uterus.

IMPORTANT: In order to maximize tamponade effect, constant pressure must be effected between the balloon wall and the tissue surface. Pressure can be achieved and maintained by applying gentle traction to the balloon shaft, before securing to the patient's leg, or attaching to a weight, not to exceed 500 grams.

INSTRUCTIONS FOR USE

IMPORTANT: Prior to transvaginal or transabdominal placement of the Bakri Postpartum Balloon, the uterus should be free of all placental fragments, and the patient should be evaluated to ensure

that there are no lacerations or trauma to the genital tract and that the source of the bleeding is not arterial.

Transvaginal Placement

1. Determine uterine volume by direct examination or ultrasound examination.
2. Insert the balloon portion of the catheter into the uterus, making certain that the entire balloon is inserted past the cervical canal and internal ostium.
3. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.

Transabdominal Placement, Post Cesarean Section

1. Determine uterine volume by direct examination (intraoperative) or ultrasound examination (postoperative).
2. From above, via access of the Cesarean incision, pass the tamponade balloon, inflation port first, through the uterus and cervix.

NOTE: Stopcock may be removed to aid placement and reattached prior to filling balloon.

3. Have an assistant pull the shaft of the balloon through the vaginal canal until the deflated balloon base comes into contact with the internal cervical ostium.
4. Close the incision per normal procedure, taking care to avoid puncturing the balloon while suturing.

NOTE: Ensure that all product components are intact and the hysterotomy is securely sutured prior to inflating the balloon. If clinically relevant, the abdomen may remain open upon inflation of the balloon to closely monitor uterine distention and confirm the hysterotomy closure.

NOTE: If clinically relevant, a B-Lynch compression suture may be used in conjunction with the Bakri Postpartum Balloon.

Balloon Inflation

With Syringe

WARNING: Always inflate the balloon with sterile liquid. Never inflate with air, carbon dioxide or any other gas.

WARNING: The maximum inflation is 500 ml. Do not overinflate the balloon. **Overinflation of the balloon may result in the balloon being displaced into the vagina.**

NOTE: To ensure that the balloon is filled to the desired volume, it is recommended that the predetermined volume of fluid be placed in a separate container, rather than relying on a syringe count to verify the amount of fluid that has been instilled into the balloon.

1. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.
2. Using the enclosed syringe, begin filling the balloon to the predetermined volume through the stopcock.
3. Once the balloon has been inflated to the predetermined volume, confirm placement via ultrasound.

NOTE: See Fig. 1 for proper placement.

4. Apply gentle traction to the balloon shaft to ensure proper contact between the balloon and tissue surface. To maintain tension, secure the balloon shaft to the patient's leg or attach to a weight, not to exceed 500 grams.

NOTE: To maximize the tamponade effect, counter pressure can be applied by packing the vaginal canal with iodine- or antibiotic-soaked vaginal gauze.

5. Connect the drainage port to a fluid collection bag to monitor hemostasis. **NOTE:** To adequately monitor hemostasis, the balloon drainage port and tubing **may** be flushed clear of clots with sterile isotonic saline.
6. Monitor the patient continuously for signs of increased bleeding and uterine cramping.

Balloon Removal

NOTE: The timing of balloon removal should be determined by the attending clinician upon evaluation of the patient once bleeding has been controlled and the patient has been stabilized. The balloon may be removed sooner upon the clinician's determination of hemostasis. The maximum indwell time is 24 hours.

1. Remove tension from the balloon shaft.
2. Remove any vaginal packing.
3. Using an appropriate syringe, aspirate the contents of the balloon until fully deflated. **The fluid may be removed incrementally to allow periodic observation of the patient.**

NOTE: In an emergent situation the catheter shaft may be cut to facilitate more rapid deflation.

4. Gently retract the balloon from the uterus and vaginal canal and discard.
5. Monitor patient for signs of bleeding.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCE

Y. N. Bakri, A. Amri, F. Abdul Jabbar: "Tamponade balloon for obstetrical bleeding," International Journal of Gynecology and Obstetrics, 74 (2001) 139-142

Vitthala S, Tsoumpou I, Anjum ZK, Aziz NA "Use of the Bakri balloon in post-partum haemorrhage: a series of 15 cases" Austr NZ J Obstet Gynaecol; August 2009 49(4):445

Gronvall M, Tikkanen M, Tallberg E, Paavonen J, Stefanovic V: "Use of Bakri balloon tamponade in the treatment of postpartum hemorrhage: A series of 50 cases from a tertiary teaching hospital" Acta Obstet Gynecol Scand. August 2012-09-18

John M. O'Brien, MD, Wendy L. Nelson, MD: "The uterine sandwich for persistent uterine atony: combining the B-Lynch compression suture and an intrauterine Bakri balloon," American Journal of Obstetrics and Gynecology, Volume 196, Issue 5 (2007) e9-e10

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LITTERATUR

- Y. N. Bakri, A. Amri, F. Abdul Jabbar: "Tamponade balloon for obstetrical bleeding," *International Journal of Gynecology and Obstetrics*, 74 (2001) 139-142
- Vitthala S, Tsoumpou I, Anjum ZK, Aziz NA "Use of the Bakri balloon in post-partum haemorrhage: a series of 15 cases" *Austr NZ J Obstet Gynaecol*; August 2009 49(4):445
- Gronvall M, Tikkanen M, Tallberg E, Paavonen J, Stefanovic V: "Use of Bakri balloon tamponade in the treatment of postpartum hemorrhage: A series of 50 cases from a tertiary teaching hospital" *Acta Obstet Gynecol Scand*. August 2012-09-18
- John M. O'Brien, MD, Wendy L. Nelson, MD: "The uterine sandwich for persistent uterine atony: combining the B-Lynch compression suture and an intrauterine Bakri balloon," *American Journal of Obstetrics and Gynecology*, Volume 196, Issue 5 (2007) e9-e10

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Vitthala S, Tsoumpou I, Anjum ZK, Aziz NA "Use of the Bakri balloon in post-partum haemorrhage: a series of 15 cases" *Austr NZ J Obstet Gynaecol; August 2009 49(4):445*

Gronvall M, Tikkanen M, Tallberg E, Paavonen J, Stefanovic V: "Use of Bakri balloon tamponade in the treatment of postpartum hemorrhage: A series of 50 cases from a tertiary teaching hospital" *Acta Obstet Gynecol Scand. August 2012-09-18*

John M. O'Brien, MD, Wendy L. Nelson, MD: "The uterine sandwich for persistent uterine atony: combining the B-Lynch compression suture and an intrauterine Bakri balloon," *American Journal of Obstetrics and Gynecology, Volume 196, Issue 5 (2007) e9-e10*

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REFERENS

Y. N. Bakri, A. Amri, F. Abdul Jabbar: "Tamponade balloon for obstetrical bleeding," *International Journal of Gynecology and Obstetrics*, 74 (2001) 139-142

Vitthala S, Tsoumpou I, Anjum ZK, Aziz NA "Use of the Bakri balloon in post-partum haemorrhage: a series of 15 cases" *Austr NZ J Obstet Gynaecol*; August 2009 49(4):445

Gronvall M, Tikkanen M, Tallberg E, Paavonen J, Stefanovic V: "Use of Bakri balloon tamponade in the treatment of postpartum hemorrhage: A series of 50 cases from a tertiary teaching hospital" *Acta Obstet Gynecol Scand*. August 2012-09-18

John M. O'Brien, MD, Wendy L. Nelson, MD: "The uterine sandwich for persistent uterine atony: combining the B-Lynch compression suture and an intrauterine Bakri balloon," *American Journal of Obstetrics and Gynecology*, Volume 196, Issue 5 (2007) e9-e10

This symbol on the label indicates that this device contains phthalates. Specific phthalates contained in the device are identified beside or below the symbol by the following acronyms:

- **BBP: Benzyl butyl phthalate**
- **DBP: Di-n-butyl phthalate**
- **DEHP: Di(2-ethylhexyl) phthalate**
- **DIDP: Diisodecyl phthalate**
- **DINP: Diisononyl phthalate**
- **DIPP: Diisopentyl phthalate**
- **DMEP: Di(methoxyethyl) phthalate**
- **DNOP: Di-n-Octyl phthalate**
- **DNPP: Di-n-pentyl phthalate**

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.