


Zilver PTX Global Data Analysis Highlights DES Benefits in Challenging Patient Populations

*Michael D. Dake, MD
University of Arizona
Tucson, AZ*

Zilver PTX Global Data Analysis




Speaker Name

Patient and Lesion Characteristics

Characteristic	Value
Patients / Lesions	2374 / 2686
Age (years)	69.5 ± 9.6
Male	71%
Diabetes	49%
Lesion length	11.1 ± 8.7 cm
Total occlusions	36%
In-stent restenosis	13%
Rutherford 4-6 (CLI)	14.7%
No patent runoff vessels	2.1%

Speaker Name

Freedom from TLR from Global Data Analysis



- High rate of freedom from TLR
- Nearly 2400 patients
- 5 year results

Speaker Name

What Patient and Lesion Factors Impact TLR Rates?

Patient Factors	Lesion Factors
Male vs. Female	Lesion length ($\leq 15\text{cm}$ vs. $>15\text{cm}$)
Diabetes	Total occlusion
Smoking Status	In-stent restenosis
Rutherford (0-4 vs. 5-6)	Calcification (none-moderate vs. severe)
Patent runoff vessels (0-1 vs. 2+)	Reference vessel diameter ($<5\text{mm}$ vs. $\geq 5\text{mm}$)
	Stent diameter
	Stent oversizing ($\leq 30\%$ vs. $>30\%$)

Speaker Name


What Patient and Lesion Factors Have the Most Impact on TLR Rates Over Time?

- Step 1:
Classification Tree analysis to rank factors for impact on TLR
 - Most impact
 - Moderate impact
 - Minimal impact
- Step 2:
Cox model to confirm impact of the patient and lesion factors and evaluate interactions

Speaker Name

Step 1: Classification Tree Analysis

- Ranks variables based on impact on TLR
 - Most impactful variable selected at each bifurcation
 - Identifies overall most and least impactful variables
- Provides input to other predictive models



Speaker Name

Step 1: Classification Tree Analysis

Impact on 1-year TLR rate

Most impact on TLR	Moderate impact on TLR	Minimal impact on TLR
In-stent restenosis	Stent diameter	Rutherford
Lesion length	Smoking status	Male vs. Female
Reference vessel diameter	Patent runoff vessels	Calcification
Total occlusion	Stent oversizing	Diabetes

Speaker Name

Step 2:

Cox Proportional Hazards Model

- Factors from Classification Tree were included in Cox Proportional Hazards model
 - All 2-, 3-, and 4-way interactions were tested
- Cox model confirmed the impact of the various patient and lesion factors
 - Identified significant interactions between some factors

Factors with Most Impact on TLR

Factors with most impact on TLR from Classification Tree	p-value ¹	Interpretation
In-stent restenosis	< 0.001	<ul style="list-style-type: none">• ISR is significant• Lesion length is significant in non-ISR lesions• Shorter lesions (≤ 15 cm) more likely to be free from TLR
Lesion length	0.34	
ISR and lesion length	0.013	
Total occlusion	< 0.001	<ul style="list-style-type: none">• Occlusion is significant• RVD is significant in non-CTO• Larger vessels (≥ 5 mm) more likely to be free from TLR
Reference vessel diameter	0.33	
Occlusion and RVD	0.037	

¹ joint p-value; not univariate

Speaker Name

Factors with Moderate Impact on TLR

Factors with moderate impact on TLR from Classification Tree	p-value ¹	Interpretation
Stent oversizing	0.043	• Oversizing (>30%) may impact TLR • No significant impact on TLR
Smoking status	0.64	
Patent runoff vessels	0.91	
Stent diameter	0.12	

¹ joint p-value; not univariate

Speaker Name

Factors with Minimal Impact on TLR

Factors with minimal impact on TLR from Classification Tree	p -value ¹	Interpretation
Rutherford	0.065	• No significant impact on TLR
Male vs. Female	0.68	
Calcification	0.98	
Diabetes	0.22	

¹ joint p -value; not univariate

Speaker Name

Conclusions from Zilver PTX Global Data Analysis

- Nearly 2400 patients from 5 studies across 16 countries
 - Including ISR, long lesions, no runoff
- High freedom from TLR sustained through 5 years
- Consistency between Classification Tree and Cox model results
- ISR and total occlusions identified as significant factors for TLR
 - Lesion length is significant in non-ISR lesions
 - RVD is significant in non-CTO
- For Zilver PTX patients, sex, calcification, and diabetes did NOT have significant impact on TLR

Speaker Name