**COOK®**

**MEDICAL**

**EchoTip® Ultrasound Needle**

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## **INTENDED USE**

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

## **NOTES**

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

## **CONTRAINDICATIONS**

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

## **POTENTIAL COMPLICATIONS**

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

## **PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew ring locked at zero centimeter mark prior to introduction, advancement or withdrawal of the device. Failure to retract needle may result in damage to endoscope.

## **SYSTEM PREPARATION**

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer slip fitting with a stopcock on side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)
2. Prepare syringe as follows:
  - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
  - b. Turn stopcock 90° to "closed" position.
  - c. Pull plunger back until it is locked into place at desired setting, creating a suction.
  - d. Set prepared syringe aside until aspiration biopsy is desired.
3. Introduce ultrasound endoscope and advance to appropriate position for biopsy.

## INSTRUCTIONS FOR USE

1. Identify desired biopsy site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance in small increments, until Luer lock fitting at base of handle meets fitting on ultrasound endoscope accessory channel. (*See fig. 2*)
3. Attach device to accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew and setting ring to desired extension, then tightening thumbscrew. **Note:** Number on operator's side of thumbscrew ring indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device has been attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to ultrasound endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing handle shaft to pre-positioned thumbscrew ring.
6. Advance needle into lesion.
7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of handle. Preserve stylet for use if additional cell collection is desired.
8. Attach Luer lock fitting of previously prepared syringe securely onto metal fitting on handle.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle, in small increments, back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection.
10. Upon completion of cell collection, retract needle into sheath by pulling back on handle shaft. Lock thumbscrew ring at zero centimeter mark.
11. Disconnect Luer lock fitting of device from accessory channel and withdraw entire device from ultrasound endoscope.
12. Unlock thumbscrew ring and extend needle. Exchange locking syringe for an empty syringe, and attach new syringe to Luer lock fitting on handle.
13. Push in plunger, then prepare aspirate specimen per institutional guidelines.
14. Aspirate remaining can be recovered for examination by flushing device.
15. For additional cell collection from same lesion, gently reinsert stylet into metal fitting on handle. **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting sheath at Luer lock fitting, advance stylet wire, in small increments, until stylet hub is engaged in fitting.

16. Additional samples may be obtained by repeating step 2 of "System Preparations" through step 14 of "Instruction for Use".

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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EchoTip is a registered trademark of Cook Urological Inc.
Cook is a registered trademark of Cook Incorporated.
**Cook Ireland Ltd.**