


MiraFlex™

MICROCATHETER

The quality you expect from
COOK embolization products,
now in a microcatheter.



COOK®

Better results start here.

We presented our engineers with a difficult challenge:

**To design a microcatheter worthy of a place in
the COOK embolization family.**



**Global
Product
Number**

**Order
Number**

G43113

MF-2.5-18-100-15

G43114

MF-2.5-18-110-15

G43115

MF-2.5-18-135-15

G43116

MF-2.5-18-150-15



Their success brings you

MiraFlex[™]
MICROCATHETER

**the only microcatheter
advanced enough to complement
COOK embolization products.**

MiraFlex embodies the latest innovations in microcatheter technology, from the same company that brought you MReye®, the world's first MR-compatible coil. Now you can deploy Nester®, Tornado® or MReye coils with the confidence that comes with a microcatheter engineered specifically for COOK coils.

It's the power to stop mixing and matching coils and catheters.
It's the power to achieve better results through better compatibility.
And it all starts with MiraFlex.

Physical attributes

- **5 durometer zones** for a distinct yet seamless transition from shaft to tip.
- **Braided construction** for added control and kink-resistance.
- **.021 inch inner diameter** permits use with a wide variety of coils, particles and spheres.
- **Lipiodol/DMSO compatible.**
- **Radiopaque band** for superior visibility.

COOK®

Diagnostic and Interventional Products

www.cookmedical.com

Distributed by:

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.

Phone: 812 339-2235, Toll Free: 800 457-4500, Toll Free Fax: 800 554-8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905 640-7110, Toll Free: 800 668-0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains

Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86