

Catheters for Peripheral Access

Cerebral

CK1
FC3
H1
H1H
H3
HN1
HN4
HN5
JB1
JB2
JB3
MAN
SIM1
SIM2
SIM3
VTK
WNBG

Visceral

C1
C2
C3
CHG2.5
CHG3
CHG-A
CHG-B
CHG-C
CHG-D
CHG-E
LEV1
MIK
RC1
RC2
RDC
RH
RIM
RUC
VS
VS1
VS2
VS3

Multipurpose

DAV
45-DEG-ANGLE
KMP
MPA
MPB
TEGT
VERT
0

Flush

NEFF
VCF
PIG
CFP
PED

Sizing

PIG-CAVA
PIG-CSC
0-CSC
PIG-WSC
0-WSC
CFP-WSC
VCF-WSC

COOK[®]

Flow Rates¹

SCBR

CATHETER		MAXIMUM FLOW RATE cc/sec	MAXIMUM RECOMMENDED PRESSURE	
French Size	Length cm	100% Oxilan® 350	kg/cm ²	psi
4.0	40	20	84	1200
4.0	65	16	84	1200
4.0	100	13	84	1200
5.0	40	18	84	1200
5.0	65	15	84	1200
5.0	70	15	84	1200
5.0	80	13	84	1200
5.0	100	11	84	1200
5.0	125	10	84	1200

HNB/HNBR

CATHETER		MAXIMUM FLOW RATE cc/sec	MAXIMUM RECOMMENDED PRESSURE	
French Size	Length cm	100% Oxilan® 350	kg/cm ²	psi
4.1	65	14	84	1200
4.1	100	10	84	1200
5.0	65	15	84	1200
5.0	100	11	84	1200

HNR

CATHETER		MAXIMUM FLOW RATE cc/sec	MAXIMUM RECOMMENDED PRESSURE	
French Size	Length cm	100% Oxilan® 350	kg/cm ²	psi
4.0	65	19	84	1200
4.0	70	20	84	1200
4.0	90	17	84	1200
4.0	100	15	84	1200
5.0	65	32	84	1200
5.0	70	32	84	1200
5.0	90	29	84	1200
5.0	100	27	84	1200
5.0	110	26	84	1200

¹Flow rates should only be used as a general guideline. Since a number of variables can produce different flow rates, please refer to product label for specific flow rate information.

COOK®

Diagnostic and Interventional Products

www.cookgroup.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812 339-2235, Toll Free: 800 457-4500, Toll Free FAX: 800 554-8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905 640-7110, Toll Free: 800 668-0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86