




## FISTULA BRUSH



Gentle preparation paves the way for successful fistula repair.

### SPECIFIC DESIGN

The only product available that is specifically designed for fistula tract identification and preparation.


### THOROUGHLY GENTLE

Cleans and gently debrides the fistula tract with a back and forth motion, removing nonvascularized tissue.

### COORDINATED SYSTEM

Pairs perfectly with the Surgisis® Biodesign™ Fistula Plug, serving as the lead to pull the plug into place.

## ORDERING INFORMATION



Used to identify, clean or debride an anal or recto-vaginal fistula tract and facilitate placement of the Surgis Biodesign Fistula Plug.

Supplied sterile in peel-open package.  
Intended for one-time use.


Not for vascular use.

Do not increase the diameter of the fistula tract.

Global Product Number	Order Number	Length (cm)
G48527	J-FB-100	46

## ACCESSORY PRODUCTS

### Surgis Biodesign Fistula Plug



Anal Fistula

Recto-Vaginal Fistula



[www.cookmedical.com](http://www.cookmedical.com)

#### COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

#### COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

#### WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

#### COOK IRELAND LTD.

O'Halloran Road, National Technology Park, Limerick, IRELAND  
Phone: 353 613 34440

© COOK 2008

SUR-BM-FB-EN-200805

AORTIC INTERVENTION

CARDIOLOGY

CRITICAL CARE

ENDOSCOPY

PERIPHERAL INTERVENTION

**SURGERY**

UROLOGY

WOMEN'S HEALTH