




# Biodesign™

## SURGISIS® POSTERIOR PELVIC FLOOR GRAFT



### Restore a pelvic floor, and a way of life.


For centuries women have suffered the discomfort and embarrassment of pelvic organ prolapse, often in silence. Today, women desire a quality of life that allows them to continue enjoying their daily activities.

Biodesign enables you to restore your patient's quality of life. Biodesign communicates with the body, signaling surrounding cells to grow across the scaffold. The result? Completely remodeled, vascularized tissue that's strong yet flexible. And unlike synthetic mesh and cross-linked biologic grafts, it is resistant to infection, encapsulation and erosion into surrounding tissue.

Biodesign is remarkably simple to handle and can easily be trimmed to any shape or size.

Cook Medical—Improving quality of life for patients through advanced tissue repair.

Global Product Number	Order Number	Size cm
<b>Posterior Pelvic Floor Graft</b>		
G48227	J-PF-POST	14 x 16
<b>4-Layer Tissue Graft</b>		
G18941	J-SLH-4S-4X7	4 x 7
G18681	J-SLH-4S-7X10	7 x 10
G18682	J-SLH-4S-7X20	7 x 20
G46612	J-PF-8X13	8 x 13
G48667	J-PF-10X15	10 x 15
G52412	J-PFV-4X7	4 x 7 (perforated)
G46604	J-PFV-9X14	9 x 14 (perforated)



Posterior Pelvic Floor Graft



4-Layer Tissue Graft

## FEATURES & BENEFITS

- Easy to hydrate and suture.
- Minimal scarring.
- Migration and erosion resistant.
- Does not encapsulate.
- Restores natural form and function.

## TIPS

- When sizing Biodesign, be sure to trim a large enough piece for the defect. Too small a piece may create tension on sutures and cause recurrence.
- Use only antiseptics and antimicrobials that have proven compatible with wound healing. Problems have been reported with gentamicin, povidone iodine, bacitracin, polymyxin B and vancomycin.
- Slow-resorbing or nonabsorbable sutures are recommended.
- After surgery, patients should abstain from heavy lifting, strenuous exercise and sexual intercourse for 4 to 6 weeks.

## DATA


1. Arunkalaivanan A. Pelvic floor reconstruction using Surgisis (porcine mesh) graft. Poster presented at: Joint Meeting of the International Continence Society and the International Urogynecological Association 34th Annual Meeting; August 25-27, 2004; Paris, France.
2. Cutner A, Spiteri M, Osborne J. Use of absorbable mesh to prevent posterior vaginal wall prolapse after laparoscopic sacrocolpopexy. *Int Urogynecol J Pelvic Floor Dysfunct.* 2001;12:577.
3. Franklin ME Jr, Treviño JM, Portillo G, et al. The use of porcine small intestinal submucosa as a prosthetic material for laparoscopic hernia repair in infected and potentially contaminated fields: long-term follow-up. *Surg Endosc.* 2008;22(9):1941-1946.
4. Hiles M, Hodde J. Tissue engineering a clinically useful extracellular matrix biomaterial. *Int Urogynecol J Pelvic Floor Dysfunct.* 2006;17(suppl 1):S39-S43.
5. Reid, RI. A retrospective analysis of recto-enterocele repair with Surgisis biomesh. Poster presented at: 31st Annual Meeting of the International Urogynecological Association; September 6-9, 2006; Athens, Greece.
6. Ross JW. Pelvic floor hernias: tissue repair or replacement? *Laparoscopy Today* Web site. 2003. [http://www.laparoscopytoday.com/2003/01/pelvic\\_floor\\_he.html](http://www.laparoscopytoday.com/2003/01/pelvic_floor_he.html). Accessed March 13, 2009.
7. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
8. Sardeli C, Axelsen SM, Bek KM. Use of porcine small intestinal submucosa in the surgical treatment of recurrent rectocele in a patient with Ehlers-Danlos syndrome type III. *Int Urogynecol J Pelvic Floor Dysfunct.* 2005;16(6):504-505.
9. Sardeli C, Bek KM, Axelsen SM. Use of porcine small intestinal submucosa mesh in the surgical treatment of recurrent pelvic organ prolapse. Poster presented at: International Continence Society 35th Annual Meeting; August 28-September 2, 2005; Montreal, Canada.
10. Sung VW, Sokol ER, Rardin CR, et al. Porcine-derived small intestinal submucosa graft-augmented rectocele repair. *J Pelvic Med Surg.* 2004;10(suppl 1):S67. Poster 124.
11. Trabuco EC, Klingele CJ, Gebhart JB. Xenograft use in reconstructive pelvic surgery: a review of the literature. *Int Urogynecol J Pelvic Floor Dysfunct.* 2007;18(5):555-563.
12. Turner RJ, Stanford EJ, Voet RL, et al. Histologic findings of porcine intestinal allograft at 2-4 years after implantation. Poster presented at: AAGL 34th Annual Meeting; November 9-12, 2005; Chicago, IL.

## PROCEDURAL AREAS

HERNIA REPAIR  
 FISTULA REPAIR  
 PLASTIC & RECONSTRUCTIVE/ENT  
 STAPLE LINE REINFORCEMENT  
 PEYRONIE'S REPAIR  
 CONTINENCE RESTORATION  
 DURAL REPAIR  
**PELVIC FLOOR REPAIR**

ANTERIOR  
 POSTERIOR  
 POSTERIOR PELVIC FLOOR GRAFT  
 4-LAYER TISSUE GRAFT

## EVOLUTION OF TISSUE REPAIR



Biodesign is a breakthrough technology that provides signals and support for the body to restore itself. It incorporates the best attributes of a biologic graft—complete remodeling and resistance to infection—yet it is easy to use, widely available and moderately priced.



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**  
 P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
 Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

**COOK (CANADA) INC.**  
 111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
 Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**  
 95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
 Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**COOK IRELAND LTD.**  
 O'Halloran Road, National Technology Park, Limerick, IRELAND  
 Phone: 353 613 34440

© COOK 2009 SUR-BM-PFLP-EN-200904

AORTIC  
INTERVENTION

CARDIOLOGY

CRITICAL  
CARE

ENDOSCOPY

INTERVENTIONAL  
RADIOLOGY

PERIPHERAL  
INTERVENTION

**SURGERY**

UROLOGY

WOMEN'S  
HEALTH