



COOK CATHETER TRAY WITH MAXIMAL STERILE BARRIERS

Order Number	Reference Part Number	Lumens	Catheter Fr	Catheter Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Fits Needle gage	Suture Needle Shape
G23246	C-UDLMY-401J-PED-IHI-CCT	2	4.0	18	5	.018	21	curved
G23247	C-UDLMY-401J-IHI-CCT	2	4.0	18	8	.018	21	curved
G23248	C-UDLMY-401J-RSC-IHI-CCT	2	4.0	18	12	.018	21	curved
G23251	C-UDLMY-501J-PED-IHI-CCT	2	5.0	16	5	.018	21	curved
G23255	C-UTLMY-501J-PED-IHI-CCT	3	5.0	16	5	.025	19	curved
G23252	C-UDLMY-501J-IHI-CCT	2	5.0	16	8	.018	21	curved
G23256	C-UTLMY-501J-IHI-CCT	3	5.0	16	8	.025	19	curved
G23257	C-UTLMY-501J-RSC-IHI-CCT	3	5.0	16	12	.025	19	curved
G23253	C-UDLMY-501J-RSC-IHI-CCT	2	5.0	16	12	.018	21	curved
G23263	C-UTLMY-701J-IHI-CCT-A-RD	3	7.0	13	15	.032	18	curved
G23265	C-UTLMY-701J-RSC-IHI-CCT-A-RD	3	7.0	13	20	.032	18	curved
G23277	C-UDLMY-801J-IHI-CCT-A	2	8.0	12	15	.032	18	curved
G23279	C-UDLMY-801J-RSC-IHI-CCT-A	2	8.0	12	20	.032	18	curved
G23284	C-UTLMY-901J-IHI-CCT-A-RD	3	9.0	11	15	.032	18	curved
G23286	C-UTLMY-901J-RSC-IHI-CCT-A-RD	3	9.0	11	20	.032	18	curved
G23299	C-UQLMY-1001J-IHI-CCT-A-RD	5	10.0	10	15	.032	18	curved
G23295	C-UQLMY-1001J-RSC-IHI-CCT-A-RD	5	10.0	10	20	.032	18	curved

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.
Additional lengths and tray offerings are available.



STANDARD COOK CATHETER TRAY

Order Number	Reference Part Number	Lumens	Catheter Fr	Catheter Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Fits Needle gage	Suture Needle Shape
G26930	C-UDLMY-401J-PED-CCT	2	4.0	18	5	.018	21	curved
G26927	C-UDLMY-401J-CCT	2	4.0	18	8	.018	21	curved
G26932	C-UDLMY-401J-RSC-CCT	2	4.0	18	12	.018	21	curved
G26936	C-UDLMY-501J-PED-CCT	2	5.0	16	5	.018	21	curved
G29553	C-UTLMY-501J-PED-CCT	2	5.0	16	5	.025	19	curved
G26934	C-UDLMY-501J-CCT	2	5.0	16	8	.018	21	curved
G26941	C-UTLMY-501J-CCT	3	5.0	16	8	.025	19	curved
G26937	C-UDLMY-501J-RSC-CCT	2	5.0	16	12	.018	21	curved
G26945	C-UTLMY-501J-RSC-CCT	3	5.0	16	12	.025	19	curved
G23259	C-UTLMY-701J-CCT-A-RD	3	7.0	13	15	.032	18	curved
G23260	C-UTLMY-701J-RSC-CCT-A-RD	3	7.0	13	20	.032	18	curved
G23273	C-UDLMY-801J-CCT-A	2	8.0	12	15	.032	18	curved
G23275	C-UDLMY-801J-RSC-CCT-A	2	8.0	12	20	.032	18	curved
G23280	C-UTLMY-901J-CCT-A-RD	3	9.0	11	15	.032	18	curved
G23281	C-UTLMY-901J-RSC-CCT-A-RD	3	9.0	11	20	.032	18	curved
G23289	C-UQLMY-1001J-CCT-A-RD	5	10.0	10	15	.032	18	curved
G23290	C-UQLMY-1001J-RSC-CCT-A-RD	5	10.0	10	20	.032	18	curved

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details. Additional lengths and tray offerings are available.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

