

We're on the path to safer connections.



Types of compatible Cook Medical enteral feeding tubes



Balloon retention gastrostomy tube



Nasoenteric feeding tube



Jejunostomy catheter



PEG tube



Locking-loop feeding tube



PEG jejunal feeding tube



Bolus PEG tube adapter

ENFit® TRANSITION CONNECTOR

We're ready when you're ready.

Cook Medical is adopting the new ISO* standards specified for enteral feeding (ISO 80369-3). These standards are also recommended by the FDA. During the transition, we are making the ENFit Transition Connector available so that

you can use our current feeding tubes with new ENFit syringes and administration sets that are already on the market. Learn more about the ENFit connector changes at stayconnected.org.

*International Organization for Standardization




ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission.



ENFit Transition Connector

- Enhances patient safety by helping to reduce harmful and potentially life-threatening misconnections between medical devices*
- Aligns with global, industry-wide standards specified by ISO and recommended by the FDA
- Ordered separately from feeding tubes, as needed
- Allows Cook feeding tubes to connect safely with administration sets and syringes that have ENFit connectors



Order Number	Reference Part Number	Quantity
G44989	FTA-ENF-30	Box of 30



ENFit®

*GEDSA (The Global Enteral Device Supplier Association) is a non-profit trade association established in part to help introduce the new ISO standard connector and facilitate adoption of the new ENFit connectors within the healthcare community. Learn more at stayconnected.org.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
 Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
 Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



CC, URO-8.5X11