



Wayne

PNEUMOTHORAX CATHETER




cookmedical.com

Wayne Pneumothorax Catheter Sets and Trays

Intended for the relief of simple, spontaneous, iatrogenic, and tension pneumothorax. The catheter is made of radiopaque polyurethane.

- The radiopaque polyurethane material is kink resistant, and the catheter features multiple sideports to maximize drainage capabilities.
- The small-bore pigtail catheter design helps reduce trauma to internal anatomy.
- The Seldinger bedside set and tray include echogenic technology for convenience, consistency, and accuracy.




Trocar method - Set

Components

- Radiopaque polyurethane pigtail catheter
- Stainless steel needle obturator
- Connecting tubing
- #11 scalpel
- One-way chest drain valve
- Three-way stopcock

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle Obturator gage	Needle Obturator Length cm	Needle Tip Configuration	Connecting Tube Length cm
G56532	C-UTPT-1020-WAYNE-IMH	10.2	29	20	19	16.5	trocar	30
G56533	C-UTPT-1400-WAYNE-IMH	14	29	20	19	16.5	trocar	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Seldinger method - Set

Components

- Radiopaque polyurethane pigtail catheter
- Pigtail straightening obturator
- Connecting tubing
- #11 scalpel
- Wire guide
- One-way chest drain valve
- 14 Fr dilator
- Three-way stopcock
- 18 gage echogenic and percutaneous entry needles
- Needle holder cup

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle Length cm	Connecting Tube Length cm	Wire Guide Diaeter inch	Wire Guide Length cm
G56535	C-UTPT-1400-WAYNE-112497-IMH	14	29	19	5*,9	30	.038	70

* This needle features echogenic technology for enhanced visualization.



Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Seldinger method - Tray

Components

- Radiopaque polyurethane pigtail catheter
- Pigtail straightening obturator
- Connecting tubing
- #11 scalpel
- 14 Fr dilator
- J-tip .038 inch diameter, 70 cm long wire guide
- 22 and 25 gage needles
- One-way chest drain valve
- 3 mL and 6 mL syringes
- Three-way stopcock
- Tinted ChloraPrep® 3 mL applicator
- Gauze sponges
- 18 gage echogenic percutaneous entry needle
- 5% lidocaine
- Fenestrated drape
- Suture with needle
- CSR wrap
- Needle holder cup

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Connecting Tube Length cm
G56537	C-UTPTY-1400-WAYNE-112497-IMH	14	29	19	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

