

Fusion[™] NEWS



Fusion[™] Titan[™] Biliary Dilation Balloon

Improving dilation outcomes

The solution to facilitating drainage caused by bile duct strictures is the Fusion[™] Titan[™] Biliary Dilation Balloon. The low profile, minimally compliant balloon design delivers the radial force needed to disrupt tight strictures and restore antegrade bile flow. The integrated catheter design, with a Nitinol stylet, provides optimal catheter integrity, which is essential in negotiating difficult anatomy.

Fusion[™] Biliary Dilation Catheter

Responding to complicated strictures

The Fusion Biliary Dilation Catheter gives you options in the management of complicated ductal strictures. The sophisticated catheter material, with a central Nitinol core, reinforces catheter stability, while the DomeTip[™] design allows for easier access through the strictured area. The Fusion biliary dilation catheter is offered in three distinct sizes – 7, 8.5 and 10 FR – with a short 2-cm taper, which facilitates entry into higher intrahepatic lesions.