

Increase the protection and accuracy of endometrial samples.





Illustration by Lisa Clark

Tao Brush™

IUMC ENDOMETRIAL SAMPLER

Reliable

The Tao Brush is a reliable diagnostic tool for endometrial sampling in an outpatient setting. The sensitivity and specificity in identifying cancer among satisfactory specimens is 100% and 96%, respectively, when atypical hyperplasia or a more severe diagnosis is made by the Tao Brush.¹

Protected

The protective sheath ensures that the Tao Brush consistently delivers a high rate of tissue suitable for diagnosis. The sheath reduces endocervical and vaginal contamination of the sample.²

Safe

The Tao Brush is easy to use, and it is well tolerated by patients. It was found to be up to one and a half times less painful than Pipelle®.³

Pipelle is a registered trademark of CooperSurgical, Inc.



www.cookmedical.com



Tao Brush™

IUMC ENDOMETRIAL SAMPLER

The brush is easy to use,¹ and its 3.5 cm length allows generous sampling of the endometrial cavity.² A plastic sheath protects the sample from endocervical and vaginal contamination.² And the brush's high diagnostic rate proves that it is a reliable sampling tool.¹ Because this brush is well tolerated by patients, it is a viable alternative to aspiration biopsy and D&C.¹

To download the instructions for use, visit www.cookmedical.com/p/tao.

Order Number	Reference Part Number	Catheter Fr	Overall Length cm	Brush Length cm
G17023	J-ES-090500	9.0	26	3.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

1. Wu HH, Casto BD, Elsheikh TM. Endometrial brush biopsy: an accurate outpatient method of detecting endometrial malignancy. *J Reprod Med.* 2003;48(1):41-45.
2. Tao LC. Direct intrauterine sampling: the IUMC Endometrial Sampler. *Diagn Cytopathol.* 1997;17(2):153-159.
3. Yang GC, Wan LS. Endometrial biopsy using the Tao brush method. A study of 50 women in a general gynecologic practice. *J Reprod Med.* 2000;45(2):109-114.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

