


Perform thyroid biopsy with ultrasound-guided visualization.



Chiba

ECHOTIP® BIOPSY NEEDLE

High-Visibility EchoTip®

Engineered to be seen under ultrasound and facilitate accurate needle placement.

Proven Design

Hundreds of microscopic dimples on the needle surface ensure a strong hyperechoic reflection and brilliant, clearly visible images.

Efficient Sampling

Ultrasound guidance reduces the number of needle passes, providing decisive nodule sampling.



www.cookmedical.com



Learn more about ultrasound-guided fine needle aspiration.

Resources

1. Robitschek J, Straub M, Wirtz E, et al. Diagnostic efficacy of surgeon-performed ultrasound-guided fine needle aspiration: a randomized controlled trial. *Otolaryngol Head Neck Surg.* 2010;142(3):306-309.
2. Seiberling KA, Dutra JC, Gunn J. Ultrasound-guided fine needle aspiration biopsy of thyroid nodules performed in the office. *Laryngoscope.* 2008;118(2):228-231.
3. Nichols K, Wright LB, Spencer T, et al. Changes in ultrasonographic echogenicity and visibility of needles with changes in angles of insonation. *J Vasc Interv Radiol.* 2003;14(12):1553-1557.
4. Culp WC, McCowan TC, Goertzen TC, et al. Relative ultrasonographic echogenicity of standard, dimpled, and polymeric-coated needles. *J Vasc Interv Radiol.* 2000;11(3):351-358.

Order Number	Reference Part Number	gage	Length cm
G27322	DCHN-23-5.0-U	23	5
G12816	DCHN-25-5.0-U	25	5

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

AORTIC INTERVENTION
 CRITICAL CARE
 ENDOSCOPY
 INTERVENTIONAL RADIOLOGY
 LEAD MANAGEMENT
 OTOLARYNGOLOGY
 PERIPHERAL INTERVENTION
 SURGERY
 UROLOGY
 WOMEN'S HEALTH
www.cookmedical.com

Customer Service Centers
EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

