

Mark the spot for radiation therapy.



EchoTip[®] Ultra

FIDUCIAL NEEDLE

Cook has led the way in EUS for over 15 years by focusing on your multidisciplinary approach. The latest result of that focus is the EchoTip Ultra Fiducial Needle, which has four gold markers preloaded in a flexible needle. Now you have a preloaded endoscopic ultrasound needle for delivering multiple markers to allow radiation oncologists deliver therapy to their patients. Cook Medical—Defining the future of EUS.



EchoTip® Ultra

FIDUCIAL NEEDLE

This device is used to implant fiducials under endoscopic ultrasound to radiographically mark soft tissue for future therapeutic procedures. Supplied sterile and is disposable - intended for single use only.

- Four pure gold fiducials preloaded in the distal tip
- Handle design facilitates marker placement
- No need to reload for up to four marker placements, saving valuable procedure time
- 22 gage needle provides a flexible delivery platform
- Tactile feel when marker is deployed

Order Number	Reference Part Number	Needle Gage	Sheath Size Fr	Fiducial Length mm	Fiducial Width mm	Number of Fiducials per Device	Minimum Accessory Channel mm
G03793	ECHO-22-F	22	5.2	5	.43	4	3.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4