

















Plastic Stents

BILIARY AND PANCREATIC

BILIARY STENTS	SIZES		STENT MATERIAL	SHAPE
	FR	LENGTH CM		
Compass BDS® Biliary Plastic Stent CBBSO	7	5, 10, 15	Ethyl vinyl acetate (EVA) copolymer with tantalum marker bands	
Cotton-Leung® Biliary Stent Set CLBS	5 7 8.5 10 11.5	5, 7, 9, 12 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18 3, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18	5, 8.5 Fr - polyethylene 7, 10, 11.5 Fr - EVA	
Cotton-Leung® Biliary Stent only CLSO	5 7 8.5 10 11.5	3, 5, 7, 9, 12, 15 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21 3, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 20 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21 3, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 20		
Cotton-Leung® Sof-Flex® Biliary Stent CLSO-SF	7 10	5, 7, 9, 12, 15 5, 7, 9, 12, 15	7, 10 Fr - polyurethane	
OASIS® Preloaded with Cotton-Leung® Biliary Stent OACL	7 8.5 10 11.5	5, 7, 9, 12, 15 5, 7, 9, 12, 15 5, 7, 9, 12, 15 5, 7, 9, 12, 15	7, 8.5 Fr - polyethylene 10, 11.5 Fr - EVA	
Zimmon® Biliary Stent Set ZEBD	5 6 7 8 10 11.5	3, 4, 5, 7, 9, 10, 12, 15 4, 7, 10, 12 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15, 18 5, 7, 9, 12, 15 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15 5, 7, 9, 12, 15	5, 6, 8, 10 Fr - polyethylene 7 Fr - EVA	
Zimmon® Biliary Stent ZSO	5 6 7 8 10	3, 4, 5, 7, 9, 10, 12, 15 4, 7, 9, 10 3, 4, 5, 6, 7, 9, 10, 12, 15, 16, 18 4, 5, 7, 9, 10, 12, 15 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15		
Cotton-Huibregtse® Biliary Stent Set CHBS	7 8.5 10 11.5	3, 5, 6, 7, 9, 10, 12, 15 5, 7, 9, 12, 15 3, 5, 7, 9, 10, 11, 12, 15 5, 7, 9, 10, 12, 15	8.5 Fr - polyethylene 7, 10, 11.5 Fr - EVA	
Cotton-Huibregtse® Biliary Stent CHBSO	7 8.5 10 11.5	4, 5, 7, 9, 10, 12, 15, 18 5, 7, 9, 10, 12, 15, 18 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 15, 18 5, 7, 8, 9, 10, 12, 15, 18		
ST-2 Soehendra Tannenbaum® Biliary Stent TTSO	8.5 10 11.5	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15	Polytetrafluoroethylene (PTFE)	
OASIS® Preloaded with ST-2 Soehendra Tannenbaum® Biliary Stent OATS	8.5 10 11.5	5, 7, 9, 12, 15 5, 7, 9, 12, 15 5, 7, 9, 12, 15		
Solus® Double Pigtail Stent with Introducer ZSS	10	1, 3, 4, 5, 7, 9, 10, 12, 15	Polyurethane with tantalum marker bands	

PLASTIC STENT SUMMARY

PANCREATIC STENTS	SIZES		STENT MATERIAL	SHAPE
	FR	LENGTH CM		
Geenen® Pancreatic Stent Set (with internal flap) GEPD	3 5 7 8.5 10 11.5	3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 3, 5, 7, 9, 12 3, 5, 7, 9, 12, 15 6, 8, 10, 12	3, 4, 5, 6, 8.5 Fr - polyethylene 7, 10, 11.5 Fr - EVA	
Geenen® Pancreatic Stent (with internal flap) GPSO	3 4 5 6 7 8.5 10 11.5	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15 2, 5, 7, 9 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 18 2, 3, 5, 7, 9, 12, 15 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 18, 20 3, 5, 7, 9, 12, 15 2, 3, 5, 7, 9, 10, 12, 15 3, 5, 7, 9, 15		
Geenen® Pancreatic Stent Set (without internal flap) GPDS	5 7	2, 3, 5, 7, 9, 12, 15 3, 5, 7, 9, 12		
Geenen® Pancreatic Stent (without internal flap) GPSOS	5 7	2, 3, 4, 5, 7, 9, 12, 15 2, 3, 5, 7, 9, 10, 12, 15		
Geenen® Sof-Flex® Pancreatic Stent (with internal flap) GPSO-SF	5	3, 5, 7, 9, 12		
Geenen® Sof-Flex® Pancreatic Stent (without internal flap) GPSOS-SF	5	3, 5, 7, 9, 12		
Zimmon® Pancreatic Stent (Radiopaque Band (RB) - Except 4 Fr) SPSOF	4 5 6 7 8.5 10	2, 4, 6, 8, 10, 12, 15, 18 2, 2.5, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 15 2, 2.5, 4, 6, 8, 10, 12, 15 2, 2.5, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 15, 18 7, 10, 12 3, 4, 5, 6, 7, 8, 10, 12, 15	4 Fr - polyethylene 5 Fr - polyethylene with Iridium marker bands 7 Fr - EVA with Tantalum marker bands 6, 8.5, 10 Fr - polyethylene with Tantalum marker band	
Zimmon® Pancreatic Stent (Radiopaque Band (RB) - Except 3 Fr) SPSOS	3 5 7	4-N, 6-N, 8-N, 10-N, 12-N, 18-N 3, 5, 7 3, 5	3 Fr - polyethylene 5 Fr - polyethylene with Iridium marker bands 7 Fr - EVA with Tantalum marker bands	
Zimmon® Pancreatic Stent Set (with ductal flap) ZEPDF	5 7	2, 4, 5, 6, 7, 8, 9, 10, 12 2, 4, 5, 6, 7, 8, 9, 10, 12, 15	5 Fr - polyethylene 7 Fr - EVA	
Zimmon® Pancreatic Stent (with ductal flap) ZPSOF	5 7	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15 2, 4, 6, 7, 8, 9, 10, 12, 15		
Zimmon® Pancreatic Stent Set (without ductal flap) ZEPDS	5 7	2, 4, 6, 8, 10, 12 4, 6, 8, 10, 12		
Zimmon® Pancreatic Stent (without ductal flap) ZPSOS	5 7	3, 5, 7, 9, 12 4, 6, 8, 10, 12, 25-NP		
Johlin Pancreatic Wedge Stent and Introducer JPWS	8.5 10	8, 10, 12, 14, 16, 18, 20, 22 8, 10, 12, 14, 16, 18, 20, 22	Polyurethane	
				

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details. RB = Radiopaque Band