


Soft on tissue,  
strong on results.



Geenen<sup>®</sup>

SOF-FLEX<sup>®</sup> PANCREATIC STENT

The Geenen Sof-Flex Pancreatic Stent, made of a strong yet flexible material, is designed to gently take the shape of the duct instead of forcing the duct to take the shape of the stent. Available with and without internal ductal flaps to meet all of your clinical needs.



[www.cookmedical.com](http://www.cookmedical.com)



## Geenen®

SOF-FLEX® PANCREATIC STENT

Used to drain obstructed pancreatic ducts.

Supplied sterile and is disposable - intended for single use only.

- Spiral side holes for enhanced drainage
- Contoured shape facilitates proper positioning
- Offset flaps minimize migration
- Sof-Flex is a "soft" pliable radiopaque material that enhances fluoroscopic visibility

Order Number	Reference Part Number	French Size Fr	Length cm	Minimum Accessory Channel mm	Wire Guide Diameter inch
<b>Internal Flaps</b>					
G49663	GPSO-SF-5-3	5	3	2.0	.035
G49664	GPSO-SF-5-5	5	5	2.0	.035
G49665	GPSO-SF-5-7	5	7	2.0	.035
G49666	GPSO-SF-5-9	5	9	2.0	.035
G49667	GPSO-SF-5-12	5	12	2.0	.035
<b>No Internal Flaps</b>					
G49668	GPSOS-SF-5-3	5	3	2.0	.035
G49669	GPSOS-SF-5-5	5	5	2.0	.035
G49670	GPSOS-SF-5-7	5	7	2.0	.035
G49671	GPSOS-SF-5-9	5	9	2.0	.035
G49672	GPSOS-SF-5-12	5	12	2.0	.035

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

