


Place biliary stents with
ease and accuracy.



Compass BDS™
BILIARY PLASTIC STENT



www.cookmedical.com



Compass BDS™

BILIARY PLASTIC STENT

This device is used to drain obstructed biliary ducts.

Supplied sterile and is disposable - intended for single use only.

- Bidirectional stent offers dual tapered ends to facilitate ease of loading and introduction into the ductal system
- Double pigtail design features 17 drainage points for each pigtail to optimize drainage
- Enhanced fluoroscopic visualization for positioning with two radiopaque bands

Order Number	Reference Part Number	French Size Fr	Length cm	Minimum Accessory Channel mm	Wire Guide Diameter inch
G55519	CBSO-7-5	7	5	2.8	.035
G55520	CBSO-7-10	7	10	2.8	.035
G55521	CBSO-7-15	7	15	2.8	.035

Recommended for use with FS-PC-7 or PC-7.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

