


More Options in Lithotripsy.



Conquest TTC[®]

LITHOTRIPTOR CABLE

We now offer the Conquest TTC in two flexible cable sizes to give you options for your lithotripsy procedures. And you can use the Conquest with a wide selection of soft wire extraction baskets as well as with our popular Soehendra Lithotripter Handle. Cook Medical: Giving you the choices you need to help you achieve quality, cost-effective patient care.



www.cookmedical.com



Conquest TTC[®]

LITHOTRIPTOR CABLE

Used to mechanically crush stones in the biliary duct when other methods of endoscopic removal have failed.

Supplied non-sterile and is disposable - intended for single use only. The "through the channel" lithotriaptor cable accepts baskets compatible with the Soehendra Lithotriaptor Handle which is reusable.

- Coil-spring cable provides flexibility and length for use through the channel of flexible endoscopes
- Flush adapter allows for contrast injection during procedure
- Compatible with Soehendra Lithotriaptor compatible soft wire extraction baskets

Order Number	Reference Part Number	French Size Fr	Length cm	Minimum Channel Size mm	Remarks
G25472	TTCL-1	8.5	170	3.2	Cable
G25518	TTCL-10	10	170	3.7	Cable
G21860	SLH-1	N/A	N/A	N/A	Handle

Soehendra Lithotriaptor Handle sold separately. Compatible with Memory 7 Fr Soft Wire, The Web[®] Extraction Baskets and Fusion[®] Lithotripsy Compatible Extraction Basket, sold separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

