

# Take complete control of your PEG placement procedures.



## Percutaneous Endoscopic Gastrostomy System

From the moment you open the Cook PEG kit, you'll notice the difference. The components are conveniently arranged in the sequential order in which they will be used, making the procedure flow as smoothly and efficiently as the nutrients you are delivering to your patients. The Cook PEG kit - available in push or pull methods - delivers the efficiency you need and the nutrients your patients need.



[www.cookmedical.com](http://www.cookmedical.com)

## Percutaneous Endoscopic Gastrostomy System

Used for percutaneous endoscopic placement to provide enteral nutrition to patients requiring nutritional support.

Supplied sterile and is disposable – intended for single use only.

- Kit components, including hemostat and tube clamp, are conveniently packaged for sequential use during the procedure for increased efficiency
- Clear feeding tube – made of flexible, high-grade medical silicone – is optimal length to give physician control throughout placement
- Sharper scissors provide for a better cut during the placement procedure
- Reference “X mark” to specify where to cut

Order Number	Reference Part Number	Feeding Tube Fr	Feeding Tube Inner Diameter mm	Insertion Wire Diameter inch	Insertion Wire Length cm
G22637	FLOW-20-PUSH-S	20	4.47	.035	260
G22638	FLOW-20-PULL-S	20	4.47	n/a	150
G22635	PEG-24-PUSH-S	24	5.5	.035	260
G22636	PEG-24-PULL-S	24	5.5	n/a	150
G31543	FLOW-20-PUSH-SSK-S*	20	4.47	.035	260
G56403	FLOW-20-PULL-SSK-S-T*	20	4.47	n/a	150
G31548	PEG-24-PUSH-SSK-S*	24	5.5	.035	260
G56402	PEG-24-PULL-SSK-S-T*	24	5.5	n/a	150
G22762	FLOW-20-PUSH-I-S**	20	4.47	.035	260
G22761	FLOW-20-PULL-I-S**	20	4.47	n/a	150
G22758	PEG-24-PUSH-I-S**	24	5.5	.035	260
G22757	PEG-24-PULL-I-S**	24	5.5	n/a	150

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

KIT CONTENTS:	
1 Fenestrated Drape	1 Silicone Feeding Tube
1 150 cm Looped Insertion Wire (Pull Only)	1 Universal Adapter
1 260 cm Wire guide (Push Only)	5 Gauze pads, 4 x 4"
1 Cold Snare	1 Scissors
	1 Hemostat
	1 Bolster
	1 Tube Clamp
	1 Bolus Adapter
	1 Twist Lock
	1 Cable Tie
	1 Bolus Feeding Adapter – Male (International Only)
	1 Patient Care Manual
	1 #11 Scalpel*
	1 Trocar Needle*
	2 Water Soluble Lubricant*
	1 25 gage x 1½" Needle**
	1 ChloroPrep® Swabs*
	1 Povidone Ointment*
	1 Xylocaine 1%, 5.0 ml glass vial**
	1 Syringe*,**

\*Safety components available within Safety Sharp Kits (-SSK) \*\*Kits without pharmaceuticals (-I Kits)

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

