


Choose the slim, smooth, in-line balloon shaft that is **easy to place** and **comfortable** for your patient.

The feeding tube sizes match perfectly with the Entuit® Initial Placement Gastrostomy Sets to streamline placement procedures.

The unique, integrated design of the in-line balloon shaft on 14 Fr and larger tubes reduces the need for overdilation during both placement and replacement procedures.






Tubes are clearly marked with fill volume and Fr size.

Silicone material is soft for patient comfort.

Balloon material is designed for durability at appropriate fill volumes.

The snug bolster stays in place.

Inner lumen diameter stays the same size even when the balloon is inflated.




Entuit® Gastrostomy BR

BALLOON RETENTION FEEDING TUBE

Entuit® Gastrostomy BR

BALLOON RETENTION FEEDING TUBE

Intended to provide gastric access for enteral feeding, medication administration, and decompression through an established gastrointestinal stoma tract. Indicated for use in percutaneous placement of an enteral feeding tube in adult and pediatric patients that require enteral feeding, medication administration, or decompression through an established gastrointestinal stoma tract.



Order Number	Reference Part Number	Fr	Marked Length cm	Balloon Volume mL	Used with Entuit Start Initial Placement Gastrostomy Set* sizes
G35614	SBRD-12	12	10	5-7	small
G35615	SBRD-14	14	10	5-7	small
G35616	SBRD-16	16	10	5-7	medium
G35617	SBRD-18	18	10	5-10	medium
G35618	SBRD-20	20	10	20	large
G35619	SBRD-22	22	10	20	large
G35620	SBRD-24	24	10	20	large

*Must be ordered separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
 Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
 Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-85X11