


Neff Percutaneous Access Set

Used for single-puncture percutaneous access to facilitate placement of a .038 inch diameter working wire guide for interventional radiology procedures.

- Sets that feature a radiopaque band help identify the location of the sheath's distal tip for accurate positioning.
- Outer sheath permits advancement of a .038 inch diameter wire guide, while the .018 inch diameter wire guide remains in place as a safety wire.
- Sets featuring hydrophilic coating reduce friction during insertion.
- The set that includes an EchoTip® needle enhances visibility under ultrasound.

Components



Order Number	Reference Part Number	Sheath ID/OD Fr	Sheath Length cm	Access Needle gage/Length cm	Needle Tip Configuration
Stainless Steel Wire Guide					
G04267	NPAS-100	4/7	18	22/15	Chiba
G04593	NPAS-104	4/7	18	21/15	Trocar
Nitinol Wire Guide					
G08426	NPAS-100-NT	4/7	18	22/15	Chiba
Nitinol Wire Guide with Radiopaque Band and Hydrophilic Coating					
G10544	NPAS-100-RH-NT	4/6	18	22/15	Chiba
G10546	NPAS-104-RH-NT	4/6	18	21/15	Trocar
G10797	NPAS-105-RH-NT-U*	4/6	18	22/20	EchoTip Chiba

*Features an EchoTip needle to enhance visibility.

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

