




Perform thyroid biopsy with ultrasound-guided visualization.



Chiba

ECHOTIP® BIOPSY NEEDLE

High-Visibility EchoTip®

The EchoTip is engineered to be seen under ultrasound and to facilitate accurate needle placement.

Proven Design

Hundreds of microscopic dimples on the needle surface ensure a strong hyperechoic reflection and brilliant, clearly visible images.

Efficient Sampling

Ultrasound guidance reduces the number of needle passes, which provides decisive nodule sampling.



www.cookmedical.com



Chiba

ECHOTIP® BIOPSY NEEDLE

To find out more about ultrasound-guided fine needle aspiration, see resources below.

- Robitschek J, Straub M, Wirtz E, et al. Diagnostic efficacy of surgeon-performed ultrasound-guided fine needle aspiration: a randomized controlled trial. *Otolaryngol Head Neck Surg.* 2010;142(3):306-309.
- Seiberling KA, Dutra JC, Gunn J. Ultrasound-guided fine needle aspiration biopsy of thyroid nodules performed in the office. *Laryngoscope.* 2008;118(2):228-231.
- Nichols K, Wright LB, Spencer T, et al. Changes in ultrasonographic echogenicity and visibility of needles with changes in angles of insonation. *J Vasc Interv Radiol.* 2003;14(12):1553-1557.
- Culp WC, McCowan TC, Goertzen TC, et al. Relative ultrasonographic echogenicity of standard, dimpled, and polymeric-coated needles. *J Vasc Interv Radiol.* 2000;11(3):351-358.

Order Number	Reference Part Number	Gage	Length cm
G27322	DCHN-23-5.0-U	23	5
G12816	DCHN-25-5.0-U	25	5

Some products or part numbers may not be available for sale in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

