

# Cook Medical Sialendoscopy Compatibility Chart

Outer Diameter (mm)	5 Fr Kolenda Introducer* G24132	6 Fr Kolenda Introducer* G24129	0.015" Wire Guide G02274	0.018" Wire Guide G02255	NCircle® Basket G23797	NGage® Basket G23798	Advance® Salivary Duct Balloon Catheter Set G34991, G34992, G34993, G34994
---------------------	------------------------------------	------------------------------------	-----------------------------	-----------------------------	---------------------------	-------------------------	---

## STORZ ERLANGEN Miniature Endoscopes

Miniature Straight Forward Telescope 0°	0.8	●	●				◐
	1.1	●	●	●			◐
	1.6	●	●	●	●	●	

## STORZ MARCHAL Miniature Endoscopes

Miniature Straight Forward Telescope 0°	0.75	●	●				◐
	0.89	●	●				◐
	1.1	●	●	●			◐
	1.3	●	●	●	●	●	◐
	1.6	●	●	●	●	●	
Examination Sheath	1.1	●	●				◐

● = Compatible with scope

◐ = Compatible alongside a scope through the 6 Fr Kolenda Introducer Sheath

\* Compatible for scope through the Kolenda Introducer Sheath