




ENTERAL FEEDING TUBE ADAPTER



Enteral feeding
through flexible design.

ADAPTER FEATURES

- Supplied with catheters used for gastroenteric feeding
- Dual adapter design ensures proper catheter configuration, eliminating inadvertent misuse
- Larger adapter facilitates feeding; smaller adapter allows flushing of catheter
- Luer lock hub provides secure connection
- Adapter design complies with AAMI standards

Ordering Information



LUER FITTING



FEEDING ADAPTER

ADAPTER USE

- Connect luer fitting for flushing of catheter.
- Connect feeding adapter for infusion of feeding solution.
- Adapters can be exchanged at any time to accommodate caregiver's needs.

Supplied sterile in peel-open package. Intended for one-time use.

Global
Product
Number

Order Number

G50006

FTA-LMLLA-R


G50007

FTA-CONV

EASY MATCH SYSTEM



Red-on-red matching system allows easy identification.



White-on-white luer fitting identifies older design.



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.

Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains

Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008

PI-BM-EFAPU-EN200803