




Experience ultra-smooth insertion with our Transitionless-Tip™.



Micropuncture®

ACCESS SET

Fewer complications

In a study by Shah et al, access-site bleeding complication rates dropped from 11.7% to 2.6% when Micropuncture was used.¹

Smooth entry

The coaxial catheter pair's Transitionless-Tip™ has been engineered to glide into the vessel with less force.²

High visibility

- The sharp, 21 gage needle features EchoTip® technology that enhances ultrasound visibility.
- The supportive nitinol wire has a palladium tip for excellent visibility under fluoroscopy.



Micropuncture®

ACCESS SET

Order Number	Reference Part Number	Introducer OD Fr	Introducer Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Needle gage/ Length cm
Standard						
G56201	MPIS-402-NT-U-SST	4.0	10	.018	40	21/7
G56202	MPIS-502-NT-U-SST	5.0	10	.018	40	21/7
Stiffened Cannula						
G35550	MPIS-402-10.0-SC-NT-U-SST	4.0	10	.018	40	21/7
G35551	MPIS-502-10.0-SC-NT-U-SST	5.0	10	.018	40	21/7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

1. Shah H, Bansal A, Srinivasan V, et al. Vascular access site bleeding complications with Micropuncture needle in patients undergoing cardiac catheterization: a comparison with conventional Cook needle. Poster presented at: The Society for Cardiovascular Angiography and Interventions 29th Annual Scientific Sessions; May 10-13, 2006; Chicago, IL. Abstract D-45.

2. Data on file.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

